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**STEP BY STEP PREPARATION OF THE IMPLANT SITE**

The actual site for the implant should be prepared in a clearly visible place for the surgeon.

View drilling sequence back page.

**Quattro implants are not recommended for bone type I.**

**HANDLING AND STERILIZATION**

Implants are supplied individually packed. Special measures are required to handle dental implants (dust avoidance and the use of gloves and titanium set of instruments).

Sterilization: The set of instruments is not sterilized and it should be before its use. The instruments should be removed from their individual plastic packaging before sterilization. During the sterilization process, strict cleanliness norms should be performed. Implants can not be resterilized.

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**IMPLANT PLACEMENT**

- Irrigate the entire implant side with sterile water or saline solution.
- Remove the vial from the blister. Then unscrew the plastic cap.
- Implant removal from the vial: For this step you must use either the Ratchet/Dynamometrical/Squared hand wrench (for manual insertion) or the c/a implant mount driver (mechanized insertion).
- Place the wrench or the c/a implant mount driver on the implant driver until it rests or fully meets the platform driver, placing partially the abutment into the prepared site. Insert the implant into the bone with a clockwise direction movement until the desired depth. Remove the driver from the implant unscrewing the driver screw using the hexagonal screwdriver of 1.25. After having removed the driver, set the cover screw to a maximum torque of 10 NW. (Cover screw\* located on its own plastic cap).
- In case of machined implant insertion it is advisable to use a limited torque.-
- Before suturing is advisable to have implant placement x-ray.
- The mucoperiosteal flap must be carefully replaced for a maximum tissue adaptation and then sutured.

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- Irrigate the entire implant site with sterile water or saline solution.
- Remove the vial from the blister. Then unscrew the plastic cap.
- Implant removal from the vial: For this step you must use the Multi-Task driver either for manual insertion with the multiconversor for c/a latch and the Ratchet/Dynamometrical/ or for mechanized insertion with c/a implant mount driver. Place the driver into the Multi-Task until the platform makes contact. Take the implant until it is partially inside the prepared site. Insert the implant into the bone with a clockwise direction movement until the desired depth. Remove the Multi-Task abutment unscrewing the screw using the hexagonal driver of 1.25. After having removed the Multi-Task abutment, set the cover screw to a maximum torque of 10 NW.
- In case of machined implant insertion it is advisable to use a limited torque.-
- Before suturing it is advisable to have implant placement x-ray.
- The mucoperiosteal flap must be carefully replaced for a maximum tissue adaptation and then sutured.
- The implant is inserted with the plastic ring. To remove the ring after implant insertion hold the post so as to leave visible the

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partially open side, Push from one of the two upper corners until the ring is removed. -In compromised situations, when the ring removal is a need, the recommended torque for the correct tightening of the ring is 5 NW/cm to avoid Morse taper effect.-

**POST OPERATIVE PROCEDURES**

The patient must be instructed to follow a routine post surgical regimen, including cold packs for 24 hours post-implantation. An antibiotic can be prescribed. Sutures can be removed after 1 week. If a removable prosthesis is used, it is highly recommended that the underside of the prosthesis be relieved and relined with a soft tissue conditioner to avoid pressure on the surgical site.

**CONTRAINDICATIONS**

Implants should not be placed if there is insufficient alveolar bone width and height to surround the bone. An insufficient bone availability and quality, poor patient oral hygiene and diseases (diabetes, etc) may interfere with osseo-integration, and result in an implant failure. You may call to our telephone number for technical assistance, get more information, or place an order. Do not use if the package is damaged.

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
V1: 19/10/2016

**LABELLING SYMBOLS**

Indications on the labels: Expiry date, Gamma-Ray sterilized, Pay attention to instructions for use, One use only, not reusable, Manufacturer, Reference, Lot, Do not use if damaged


The product remains sterile for 5 years from the time of sterilization. The expiry date of the implant is indicated on the label.

**PRODUCT COMPOSITION**


Mozo-grau manufactures dental implants and abutments made of titanium grade IV and grade V. Biocompatible materials.

V1:19/10/2016

**USO DEL VIAL / VIAL HANDLING**

**USO DEL VIAL MTA / MTA VIAL HANDLING**

**INSTRUCTIONS FOR USE**

**TREATMENT PLANNING PRIOR TO SURGERY**

Radiographs should be accurately measured to ensure the correct selection of implant's length and to avoid space between the jawbone and the sinus, the nasal cavity, the mandibular canal, or drilling the internal aspects of the mandible. There should be at least 2mm of bone under the actual implant. For example, to place a 8mm-length implant, 10mm of bone height is required.

The crest profile should be carefully examined to calculate the insertion angle to achieve parallelism with additional implants and natural dentition.

IMPORTED AND MANUFACTURED BY MOZOGRAU, S.A. (Spain)

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