INSTRUCTIONS FOR USE

PRODUCT IDENTIFICATION

Ticare dental implants are root-form, self-tapping endosseous dental implants to be surgically implanted in the jaw bone to hold and retain a dental prosthesis. There are two options with regard to prosthetic restorations: -In-Hex for internal connection -Osseous for external connection

The following dental implant variants are available:

- In-Hex
- o Internal connection
- o Internal connection Quattro
- o Internal connection Multitask (MTA)
- Osseous
- o External connection
- o External connection Quattro
- o External connection Multitask (MTA)

INDICATION FOR USE

Ticare dental implants are intended to be implanted in the maxilla or mandible jaw bone to serve as a union between the jaw bone and a dental prosthesis for partial or total replacement of teeth in edentulous patients.

Ticare dental implants are indicated for single-stage or two-stage procedures to support screw-retained restorations and can be used for immediate loading when good primary stability is achieved and with appropriate occlusal loading. Osseous Quattro and In-Hex Quattro dental implants are indicated to support permanently fixed restorations and, due to their thread design, are not indicated for use in bone type I. Small diameter (3.3mm) Ticare Mini implants are intended to replace a lateral incisor in the maxilla and/or a central or lateral incisor in the mandible and should not be used in the molar region.

CONTRAINDICATIONS

Implants should not be placed if there is insufficient alveolar bone width and height to surround the implant. Insufficient bone quantity and/or quality, poor patient oral hygiene and diseases (diabetes, etc) may interfere with osseointegration and result in implant failure. Dental implant procedures are contraindicated for patients who are medically

unfit to undergo an oral surgical procedure. Do not use if the patient presents known allergy or hypersensitivity to any of the material components – See section on Material Composition. Small diameter (3.3mm) Ticare Mini implants are not intended to be used to replace a central incisor, a canine, a premolar or a molar in the maxilla and are not intended to be used to replace a canine, a premolar or a molar in the mandible. Due to their specific thread design, Osseous Quattro and In-Hex Quattro dental implants are not indicated for use in bone type I.

GENERAL WARNINGS

Do not use if the individual package is open or damaged.

GENERAL PRECAUTIONS

Caution: Federal (US) law restricts this device to sale by or on the order of a licensed dentist or physician. Ticare dental implants are intended to be used only by qualified specialist professionals (dentist, doctors experts in dentistry or stomatology, maxillofacial surgeons).

Ticare dental implants are intended to be used only in conjunction with Ticare abutments and associated instruments according to the instructions and recommendations of Mozo Grau. Combining implants, instruments and prosthetic components that are not dimensioned for correct mating can lead to mechanical failure, damage to tissue or unsatisfactory aesthetic results.

Failure to observe the indications for use and directions for use may result in implant failure. Bone loss and/or bone resorption may lead to implant failure to osseointegrate and mechanical failure including implant fatigue fracture can occur.

TREATMENT PLANNING PRIOR TO SURGERY

Careful clinical and radiological examination of the patient is essential to treatment planning prior to surgery.

Precaution: Special care is advised when treating patients who have local or systemic risk factors that may interfere with the healing process of bone or soft tissue. These risk factors include smoking, poor oral hygiene, uncontrolled diabetes, oral radiotherapy, infections in neighbouring bone and certain therapies including steroids and chemotherapy. Bruxism or other jaw disorder may impede the success of dental implant

procedures and may require the consideration of alternative treatment options. Advise patients who smoke to stop smoking before undergoing dental implant procedures.

Precaution: Routine use of dental implants in paediatric patients is not recommended until jaw bone growth is completed. Radiographs should be accurately measured to ensure the correct selection of implant length and to avoid the space between the jawbone and sinus, the nasal cavity, the mandibular canal, or drilling the internal aspects of the mandible.

Warning: Failure to use radiographic measurements when planning surgery and the drilling sequence can result in permanent damage to vital oral structures. Drilling beyond the depth intended for lower jaw surgery may result in numbness to the lower lip or chin or lead to hemorrhage in the floor of the mouth.

There should be at least 2mm of bone under the actual implant. For example, to place an 8mm-length implant, 10 mm of bone height is required.

The crest profile should be carefully examined to calculate the insertion angle to achieve parallelism with additional implants and natural dentition.

Warning: Failure to refer to anatomic knowledge and specific pre-operative radiographs may result in damage to vital oral structures when drilling in the jaw bone.

HANDLING AND STERILIZATION

Implants are supplied sterile and individually packed. Store in a dry place away from direct sunlight at room temperature. Implant materials are stable at under normal atmospheric pressure and temperature conditions. However, to prevent damage to packaging, avoid exposing the product to temperatures above 45°C and relative humidity exceeding 70%.

Warning: Do not use if the individual package is open or damaged.

Warning: Sterile, single-use device. Do not re-use. Do not re-process. Re-use could cause cross contamination. Also, due to mechanical stresses inherent to implant placement, re-use of the dental implant may lead to implant fracture or failure. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Use of a non-sterile device may lead to infection. Use special precaution when handling dental implants. Avoid

ambient dust or contamination. Use gloves. Use titanium instruments when holding or gripping implants.

Precaution: Use aseptic technique to handle the implant once removed from its individual packaging. All instruments used should be sterilized before use. Use gloves to handle the product. Avoid dust and contamination of the surgical site.

DIRECTIONS FOR USE – OSSEOUS AND INHEX IMPLANT PLACEMENT

Prepare the implant site in a place that is clearly visible for the surgeon.

Precaution: To avoid damage to bone, bone loss and/or bone necrosis, follow the recommended surgical drilling sequence using constant irrigation and recommended drill speeds.

Precaution: Dental implants and abutments are small in size. Care must be taken to ensure that implants are not swallowed or aspirated by the patient.

1. Irrigate the entire implant site with sterile water or saline solution.

2. Remove the vial from the blister. Then unscrew the plastic cap.

3. To remove the implant from the vial, use either the Ratchet/Dynamometrical/Squared hand wrench (for manual insertion) or the contra-angle implant mount driver (for mechanized insertion).

Place the wrench or the contra-angle implant mount driver on the implant driver until it rests or fully meets the platform driver, partially placing the abutment into the prepared site. Insert the implant into the bone with a clockwise direction movement until the desired depth is reached. For mechanized implant insertion, use limited torque.

Precaution: Avoid overtightening the implant. Do not exceed the recommended maximum torque. Overtightening the implant may lead to implant damage or fracture or to bone necrosis at the bone site.

Remove the driver from the implant unscrewing the driver screw using the hexagonal 1.25 mm screwdriver.

After removing the driver, set the cover screw to a maximum torque of 10 N/cm. (Cover screw* located in its own plastic cap).

4. Before suturing, it is advisable to obtain an implant-placement x-ray.

5. Carefully replace the mucoperiosteal flap for maximum tissue adaptation and suture.

DIRECTIONS FOR USE – MTA OSSEOUS AND MTA INHEX IMPLANT PLACEMENT

Prepare the implant site in a place that is clearly visible for the surgeon.

Precaution: To avoid damage to bone, bone loss and/or bone necrosis, follow the recommended surgical drilling sequence using constant irrigation and recommended drill speeds.

Precaution: Dental implants and abutments are small in size. Care must be taken to ensure that implants are not swallowed or aspirated by the patient.

1. Irrigate the entire implant site with sterile water or saline solution.

2. Remove the vial from the blister. Then unscrew the plastic cap.

3. To remove the implant from the vial, use either the Ratchet/Dynamometrical/Squared hand wrench (for manual insertion) or the contra-angle implant mount driver (for mechanized insertion).

Place the wrench or the contra-angle implant mount driver on the implant driver until it rests or fully meets the platform driver, partially placing the abutment into the prepared site. Insert the implant into the bone with a clockwise direction movement until the desired depth is reached. For mechanized implant insertion, use limited torque.

Precaution: Avoid overtightening the implant. Do not exceed the recommended maximum torque. Overtightening the implant may lead to implant damage or fracture or to bone necrosis at the bone site.

Remove the Multi-task MTA abutment from the implant unscrewing the screw using the hexagonal 1.25 mm screwdriver.

After removing the Multi-task MTA abutment, set the cover screw to a maximum torque of 10 Nm. (Cover screw* located in its own plastic cap).

When using the MTA abutment, the implant is inserted with the plastic ring. To remove the ring after implant insertion, hold the post so that the partially open side is visible. Push from one of the upper 2 corners to remove the ring. In compromised

situations where the ring must be removed, the recommended torque for correct tightening of the ring is 5 N/cm to avoid Morse taper effect.

4. Before suturing, it is advisable to obtain an implant-placement x-ray.

5. Carefully replace the mucoperiosteal flap for maximum tissue adaptation and suture.

SIDE EFFECTS AND POST OPERATIVE PROCEDURES

Precaution: Surgical complications such as infection, mucositis and peri-implantitis as well as side-effects such as pain and swelling may occur following completion of the implant procedure.

Instruct the patient to follow a routine post-surgical regime, including cold packs for 24 hours post-implantation.

Precaution: Avoid infection in the post-operative phase. Advise the patient to maintain strict oral hygiene. Smoking increases the risk of infection. Nicotine may interfere with the healing process. Advise the patient to stop smoking when undergoing dental implant procedures.

An antibiotic can be prescribed. Sutures can be removed after one week. If a removable prosthesis is used, it is highly recommended that the underside of the prosthesis be relieved and relined with a soft tissue conditioner to avoid pressure on the surgical site.

After implant placement, careful evaluation of bone quality and primary stability will determine when implants may be loaded.

Precaution: Avoid risks of prosthetic overload and use particular caution when placing narrow implants in the posterior region.

Precaution: Lack of adequate quantity or quality of remaining bone, infection and disease may cause failure of the implant to osseointegrate both immediately after surgery or over the long-term.

Precaution: Regular patient follow-up is needed to secure long-term treatment outcome and to inform the patient about appropriate oral hygiene.

MATERIAL COMPOSITION

Ticare dental implants and dental implant abutments are made of Grade IV Titanium and Grade V Titanium Alloy.

LABELLING SYMBOLS

The following symbols are used on the labelling to convey essential information to the user:

See instructions for use

Manufacturer

Reference

Lot



CE 0051 CE Certificate

MANUFACTURER: MOZO GRAU, S.A.







C/ Santiago López González, 7
47197 VALLADOLID (SPAIN)
Tel. +34 983 211 312 fax: +34 983 304 021
e-mail: sales@ticareimplants.com
www.ticareimplants.com

Rev. 16 | 09/10/2020
Ref: APEIV20A

*** DRILLING SEQUENCE**

Profile drills should not be used for implants 6mm long

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MOZGRAU, S.A.