MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING 
	2.2.2 INSTRUCTIONS FOR USE (CASTABLES)	Rev.: 01
		Date: 15/10/2024
		Página 1 de 6

ESPAÑOL

ALCANCE

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ENGLISH

SCOPE

This instruction for use has its scope in TICARE castables manufactured by Mozo Grau.

TICARE Castable abutments are pieces made out of meltable plastic that are the support for the generation of the definitive prosthesis in the corresponding material.

A castable abutment serves as an intermediary component in dental implant procedures. The castable abutment is a prefabricated component that can be used to create a custom abutment for either a cement-retained or screw-retained prosthesis. Professionals create a custom abutment by waxing a plastic burnout pattern based on the patient's specific needs. This pattern is then cast using a lost-wax technique to form the final abutment. After implant placement, a temporary filling material covers the screw access channel, allowing easy future access and adjustments if needed.

INTENDED USE

TICARE Castable abutments are intended to be screwed onto the laboratory replica to mould on itself the shape and volume of the definitive prosthesis using dental materials.

INDICATIONS FOR USE

TICARE Castable Abutment is indicated for the following medical conditions:

- Oral endosseous implantation in the upper and lower jaw to restore chewing function.
- Aesthetic dental rehabilitation.

The castable abutment is intended for dental laboratory use only.

The castable abutments are indicated for the realization of screwed prostheses.

Cover the external surfaces of the plastic coping completely with wax (minimum wax thickness of 0.3 mm).

Correct transfer from model to patient: To ensure a correct transfer of the abutment position from the master cast to the patient, an individual index can be fabricated on the cast. In case of single crowns, the index is secured with the connection of the implant/abutment and ideally supported from the adjacent teeth as well; and in case of bridges, the abutments are splinted to one another.

To remove the prosthetic devices from the analogue or implant always use the corresponding screwdriver.

Castable UCLA abutments:

These castable abutments are pieces made out of meltable plastic that are the support for the generation of the definitive prosthesis in the corresponding material which are connected directly to the implant/analogue.

Castable posts:

These castable abutments are pieces made out of meltable plastic that are the support for the generation of the definitive prosthesis in the corresponding material which are connected to the implant/analogue through a tapered abutment.

USER

TICARE Castable abutments should only be used by trained and specialised dental and dental professionals

INTENDED TARGET GROUPS

TICARE Castable abutment is to be used in a replica of an edentulous or partially dentate adult patient's mandible.

CONTRAINDICATIONS AND WARNINGS

The castable abutment is for single-use only.

Reuse of the castable abutment can result in loss of functionality.

TICARE castable abutments are not sold sterile.

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING 
	2.2.2 INSTRUCTIONS FOR USE (CASTABLES)	Rev.: 01
		Date: 15/10/2024
		Página 5 de 6

The products must be used only for the function for which they were specified. It is essential to use the appropriate products for each type of implant, for this consult the product catalogue.

The castable abutment or screw must never be changed or modified.

Do not use if the individual package is open or damaged.

For the safe and effective use of the TICARE Castable abutments, specific training and study of appropriate textbooks is strongly recommended. THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR PROPER TRAINING. The responsibility for proper patient choice, adequate training, experience in implant placement, and provision of appropriate information for informed consent rests with the practitioner. Improper technique may result in failure and/or loss of supporting bone.

PACKAGING AND STORAGE

TICARE castable abutments do not need special conditions for storage, because the raw material is stable under atmospheric conditions of pressure and temperature. Transport and Storage temperature <55°C and humidity <75%. TICARE castable abutments are stable at atmospheric conditions of transport and storage.

The product is delivered clean, but not sterile. All information about the product is on the label. Product packaging is done individually.

Store the set at a moderate temperature in a dry, dust-free and ventilated place where corrosive vapors are not produced.

DISPOSAL

The product should be disposed of in an environmentally friendly manner in accordance with local regulations. Hazardous waste of contaminated devices or sharp objects should be disposed of in suitable containers which meet the specific technical requirements.

DETECTION AND NOTIFICATION OF INCIDENTS WITH HEALTH PRODUCTS

Important notice: Any patient or healthcare personnel who detects any serious incident related to TICARE Castables must notify the manufacturer, Mozo Grau S.A. and the competent authority of the Member State in which the user and/or patient are established.

LABELLING SYMBOLS

SYMBOLS	DESCRIPTION
	Reference of the producto
	Batch code
	Manufacture date
	Manufacturer
	Medical device
	Non sterile
	CE Mark
	Consult electronic instructions for use

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING 
	2.2.2 INSTRUCTIONS FOR USE (CASTABLES)	Rev.: 01
		Date: 15/10/2024
		Página 6 de 6

	Unique device identifier
---	--------------------------

MANUFACTURE
Mozo Grau S.A.
C/ Santiago López González, N°7
47197 Valladolid (Spain)
Tel: +34 983 309 602
e-mail: info@ticareimplants.com
www.ticareimplants.com