MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING 
	2.2.3 INSTRUCTIONS FOR USE (IMPRESSION DEVICES)	Rev.: 01
		Date: 15/10/2024
		Página 1 de 8

ESPAÑOL

ALCANCE

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ENGLISH

SCOPE

This instruction for use has its scope in TICARE impression devices manufactured by Mozo Grau.

Any prosthetic rehabilitation on implants starts with taking correct impressions, enabling precise recording of the implant location, orientation of the implant and soft tissue profile. An impression post is an accessory that records the position of the dental implant or its abutment with respect to adjacent structures; they are indexed to ensure that their location can be reproduced three-dimensionally. Impression posts (in the case of analogical impressions) or scanbodies (in the case of digital impressions) are screwed to the implants either by a screw or by its own thread. Impressions can be taken through analogue procedures using additive silicones or by digital procedures using an intraoral scanner.

There are two types of analogue impressions: Open-tray impressions, where the impression post remains inside the silicone when the tray is removed from the patient's mouth; and Closed-tray impressions, in which the impression post remains inside the patient's mouth when the tray is removed. The impression post must be removed from the mouth and manually repositioned in the impression.

For digital impressions, the scanbodies of the corresponding platform are placed on the implant with the help of the screwdriver and obtain one or several scans of the patient's mouth, with at least one of them including the Mozo Grau scanbodies. Then, a .stl archive is exported and the model is generated with the help of its specific library.

INTENDED USE

TICARE Impression devices are intended to be screwed either onto implant in mouth or onto the laboratory replica on a laboratory working model to retrieve implant 3D position and orientation and transfer it to digital software or laboratory replica.

INDICATIONS FOR USE

TICARE Impression devices are intended to be used by qualified dental surgeons or dental healthcare professionals who have received additional training in dental implant procedures. They are intended to be used in edentulous and partially dental patients, with established bone and who have finished its growth period.

TICARE Impression devices are indicated for the following medical conditions:

- Oral endosseous implantation in the upper and lower jaw to restore chewing function.
- Aesthetic dental rehabilitation.

To remove the prosthetic devices from the implant always use the corresponding screwdriver.

Regarding soft tissues, the choice of the appropriate impression coping is determined depending on the abutment chosen for the soft tissues healing and modelling. The impression coping chosen must exceed the healing abutment height and respect the emergence profile of the healing abutment used in the previous stage.

The TICARE digital flow libraries are specially designed to fit the TICARE digital impression devices.

The 3D model digital impression devices are available on the Osseous and InHex mini, standard and maxi platforms and are compatible with all the prosthetics in the TICARE catalogue.

Impression posts:

Accessories that record the position of the dental implant or its abutment with respect to adjacent structures; they are indexed to ensure that their location can be reproduced three-dimensionally by analogue techniques.

Scanbodies:

Accessories that record the position of the dental implant or its abutment with respect to adjacent structures; they are indexed to ensure that their location can be reproduced three-dimensionally by digital techniques. These devices require of a specific library. There are also scanbodies compatible with other manufacturers' devices.

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING 
	2.2.3 INSTRUCTIONS FOR USE (IMPRESSION DEVICES)	Rev.: 01
		Date: 15/10/2024
		Página 6 de 8

USER

TICARE Impression devices should only be used by trained and specialized dental and dental technicians.

INTENDED TARGET GROUPS

TICARE Impression devices is to be used in edentulous or partially dentate adult patients.

CONTRAINDICATIONS AND WARNINGS

TICARE impression devices can be reused following manufacturer instructions of cleaning/ preservation and respecting the number of uses.

TICARE impression devices should never be used if damaged or worn.

TICARE impression devices are not sold sterile. TICARE impression devices must be sterilized before use and kept within the sterile field, ensuring aseptic conditions during implant procedures.

The products must be used only for the function for which they were specified. It is essential to use the appropriate products for each type of implant or technique of impression, for this consult the product catalogue.

The impression devices or screws must never be changed or modified.

The impression devices or screws must only be used and handled by dental professionals.

The torque applied to a screw for tightening an impression device placed into the implant must not exceed the recommended torque defined by the original component manufacturer.

The torque applied to the screw must not exceed 3 Ncm.

Do not use if the individual package is open or damaged.

For proper fit and positioning of TICARE digital impression devices use only original TICARE libraries.

For the safe and effective use of the TICARE Impression devices, specific training and study of appropriate textbooks is strongly recommended. THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR PROPER TRAINING. The responsibility for proper patient choice, adequate training, experience in implant placement, and provision of appropriate information for informed consent rests with the practitioner. Improper technique may result in failure and/or loss of supporting bone.

CLEANING AND STERILIZATION

It is recommendable to clean, and sterilize TICARE impression devices before use after delivery since products are supplied non-sterile (they must be cleaned and disinfected after having removed packaging).

A proper cleaning and disinfection are essential prerequisites for proper sterilization.

Recommended cleaning method:

- Immerse devices in cleaning solution. Mozo Grau does not recommended any specific cleaning and/or disinfection agents. Only specifically formulated cleaning agents should be used. Solution with a pH value of 7-10 and a temperature of approx. 35°C/95°C.
- Rinse and flush with purified or sterile water. Mineral residues from hard water can result in staining of the device.
- Dry with compressed air or wipes.

Please do not forget that the disinfectant used in the pre-treatment is only for personal protection and does not replace the later disinfection after cleaning.

Method of sterilization:

- Place product inside sterilization bags.
- Place the product inside the autoclave. Do not sterilize different types of materials in contact.
- Sterilization in autoclave by moisty heat 134 ° C (273 ° F) - 5 minutes. Make sure elements inside the autoclave are not rusty. Only for the United States: Validated parameters are Temperature 270°F / 132°C for 4 minutes and a drying time of 20 minutes.
- The exposure temperature for products and sterilization trays shall not exceed 137 ° C.
- When the autoclave is finished, remove the product form the autoclave.

REUSABILITY INFORMATION

TICARE impression devices are reusable devices which shall be inspected before each reuse to ensure

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING 
	2.2.3 INSTRUCTIONS FOR USE (IMPRESSION DEVICES)	Rev.: 01
		Date: 15/10/2024
		Página 7 de 8

that the integrity and performance continue to be maintained. Impression copings shall be discarded if any of the following criteria are met:

- If any wear, abrasion, deformations or corrosion is visible on the component.
- If the impression coping / scanbody does not seat accurately or has a loose fit in the implant, the base, or the respective replica.
- If with light pressure the screwdriver does not engage or slips in the receptacle of the screw.

PACKAGING AND STORAGE

TICARE castable abutments do not need special conditions for storage, because the raw material is stable under atmospheric conditions of pressure and temperature. Transport and Storage temperature <55°C and humidity <75%. TICARE castable abutments are stable at atmospheric conditions of transport and storage.

The product is delivered clean, but not sterile. All information about the product is on the label. Product packaging is done individually.

No special storage conditions are needed since the raw materials used are stable under normal pressure and temperature conditions.

Store the set at a moderate temperature in a dry, dust-free and ventilated place where corrosive vapors are not produced.

DISPOSAL

The product should be disposed of in an environmentally friendly manner in accordance with local regulations. Hazardous waste of contaminated devices or sharp objects should be disposed of in suitable containers which meet the specific technical requirements.

DETECTION AND NOTIFICATION OF INCIDENTS WITH HEALTH PRODUCTS

Important notice: Any patient or healthcare personnel who detects any serious incident related to the TICARE impression devices must notify the manufacturer, Mozo Grau S.A. and the competent authority of the Member State in which the user and/or patient are established.

LABELLING SYMBOLS

SYMBOLS	DESCRIPTION
	Reference of the product
	Batch code
	Manufacture date
	Manufacturer
	Medical device
	Non sterile
	CE Mark

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING CE
	2.2.3 INSTRUCTIONS FOR USE (IMPRESSION DEVICES)	Rev.: 01
		Date: 15/10/2024
		Página 8 de 8

	Consult electronic instructions for use
	Unique device identifier

MANUFACTURE
Mozo Grau S.A.
C/ Santiago López González, N°7
47197 Valladolid (Spain)
Tel: +34 983 309 602
e-mail: info@ticareimplants.com
www.ticareimplants.com