MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING CE 0051
	2.2.2 IFU (EXPANDER/ EXPANSORES)	Rev.: 03
		Date: 05/12/2024
		Página 1 de 13

ESPAÑOL

ALCANCE

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ENGLISH

SCOPE

This instruction for use has its scope in Dental Expander manufactured by Mozo Grau.

The Dental Expander of TICARE implant systems are mechanical cutting tools, made of stainless steel. They are straight in shape and are used in conjunction with a handpiece connected to a motor or manually connected to a ratchet wrench.

All rotatory and oscillating instruments manufactured by Mozo Grau are according to the current ISO 3823 standard


INDICATIONS FOR USE

These dental instruments are indicated for the drilling and / or preparation of the bone in oral implantology processes for the placement of a dental implant.

Function of use of the Expanders: **Bone expansion or compaction system, measured and progressive**

Characteristics:


- They are stainless Steel
- Indicated for thin ridges or to compact the bone in type IV bone ridges.
- It is advisable to use the expanders sequentially, starting with number 1 and ending with the expander corresponding to the length and diameter of the selected implant.
- It is also appropriate for performing the cortical Split technique.
- Advisable for use only in the upper jaw.
- Its design and marking allows for progressive depth and promote maximum control of the depth and insertion axis.
- If excessive resistance is found, it is recommended to stop the procedure for a few seconds to avoid over compression of the bed.
- Designed with depth markings:

- They are supplied together with a template with the expansion sequence for each diameter and length of the implant. The expander template specifies, for each expander, the diameter of the socket generated at each depth mark. The template follows a colour code

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING CE 0051
	2.2.2 IFU (EXPANDER/ EXPANSORES)	Rev.: 03
		Date: 05/12/2024
		Página 9 de 13

according to the colour code of the different Ticare implant platforms: Green for 3.3mm, 3.4mm implants, blue for 3.75mm implants, purple for 4.25mm implants, turquoise for 5mm implants.

	EXPANDER Nº1	EXPANDER Nº2	EXPANDER Nº3	EXPANDER Nº4	EXPANDER Nº5	EXPANDER Nº6
	→ 8 mm: Ø2.0 mm → TIP: Ø1.0 mm	→ 15 mm: Ø3.0 mm → 13 mm: Ø2.8 mm → 11.5 mm: Ø2.65 mm → 10 mm: Ø2.45 mm → 9 mm: Ø2.3 mm → 8 mm: Ø2.2 mm → TIP: Ø1.4 mm <i>For 3.30 and 3.40 Implants</i>	→ 15 mm: Ø3.5 mm → 13 mm: Ø3.3 mm → 11.5 mm: Ø3.1 mm → 10 mm: Ø2.9 mm → 9 mm: Ø2.75 mm → 8 mm: Ø2.6 mm → TIP: Ø1.7 mm <i>For 3.75 Implants</i>	→ 15 mm: Ø4.1 mm → 13 mm: Ø3.9 mm → 11.5 mm: Ø3.6 mm → 10 mm: Ø3.4 mm → 9 mm: Ø3.25 mm → 8 mm: Ø3.1 mm → TIP: Ø2.1 mm <i>For 4.25 Implants</i>	→ 15 mm: Ø4.5 mm → 13 mm: Ø4.2 mm → 11.5 mm: Ø4.0 mm → 10 mm: Ø3.8 mm → 9 mm: Ø3.6 mm → 8 mm: Ø3.4 mm → TIP: Ø2.4 mm 	→ 15 mm: Ø4.9 mm → 13 mm: Ø4.6 mm → 11.5 mm: Ø4.3 mm → 10 mm: Ø4.1 mm → 9 mm: Ø3.9 mm → 8 mm: Ø3.8 mm → TIP: Ø2.7 mm <i>For 5.0 Implants</i>

USER

Expanders are designed to be used by qualified dental surgeons or dental health professionals who have received additional training in dental implant procedures (dental implantology).

INTENDED TARGET GROUPS

Expanders are intended to be used in edentulous and partially dentate patients, with consolidated bone and completed growth period, and who are subject to dental implant treatments or aesthetic oral rehabilitation.

CONTRAINDICATIONS

Contraindications include:

- 1. Patients allergic or hypersensitive to the raw materials used. Stainless steel
- 2. Patients with mental problems that prevent the patient's cooperation.
- 3. Patients who abuse drugs and alcohol.
- 4. Products that have changed their properties with use, such that they are not safe for the functions indicated, will not be used.

WARNINGS

The instruments used exclusively by personnel trained and specialized in dentistry. For the safety and effectiveness of the use of the system, it is strongly suggested that specific training be carried out, as well as the study of appropriate books. THESE INSTRUCTIONS ARE NOT INTENDED TO BE A SUBSTITUTE FOR ADEQUATE TRAINING.

Instrument handling must be carried out carefully, avoiding the use of damaged products as this may compromise their effectiveness. Check the condition of the cutting edges. If there is any sign of wear, discard the bone taps.

After each use, clean and disinfect the product. For good maintenance of the instruments, the cleaning and sterilization instructions must be followed. Do not store wet instruments after the cleaning and/or sterilization process.

Each instrument must be used for its intended use.

Allergen information:

The materials used for the manufacture of dental surgical burs:

- 465 KL (Stainless Steel) (C, Si, Mn, P, S, Cr, Mo, Ni, Ti, Fe)

PRECAUTIONS

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING  0051
	2.2.2 IFU (EXPANDER/ EXPANSORES)	Rev.: 03
		Date: 05/12/2024
		Página 10 de 13

All product is supplied NON-STERILE, and must be used under aseptic conditions.

The products included in this instruction can be reused. The reuse of these products is possible, as long as strict control is maintained over their cleaning and sterilization.

It is recommended that after 20 uses and their corresponding sterilizations after each surgery, the uncoated stainless steel instruments be replaced to guarantee the cutting edge.

Stainless steel is a resistant material, although it is not rust-free if it is not properly maintained.

CLEANING AND STERILIZATION

A) Treatment instructions/warnings

To avoid damage, do not use metal brushes or cleaning sponges.

Only use cleaning and disinfecting solutions with a pH value between 4.5 and 10. These instruments are supplied non-sterile and must be cleaned and sterilized before use.

B) Post processing restriction

The end of the product's useful life is normally marked by wear and tear resulting from incorrect use and handling.

C) Cleaning, maintenance and sterilization of instruments

Surgical instruments can deteriorate if proper maintenance is not carried out. At Mozo Grau we recommend our customers carry out appropriate cleaning and maintenance procedures before and after use to ensure a good useful life of the instrument.

C.1) During the use of surgical material: dirty products must be separated so as not to contaminate the instrument tray that is being used. Dirty products should be cleaned and replaced in the instrument tray before proceeding to sterilize it with the necessary components.

C.2) Treatment after surgical use:

Major dirt should be removed from the products immediately after use. To do this, use pure water or a solution with disinfectant, which should not contain aldehyde (otherwise the blood would settle into the product).

1 Selection of the product for cleaning/disinfection. It would be recommended to use a neutral/enzymatic cleaner (pH allowed between 6 and 10): It should not contain bleach, organic solvents (Alcohols, Ethers, etc.), oxidants (e.g., hydrogen peroxide), or halogens (chlorine, bromine, iodine). Be careful with tap water that is too chlorinated. Avoid solution with aldehyde, due to its ability to bind blood to the surgical material. No specific cleaning or disinfecting agent is recommended, but it is necessary to observe the instructions on concentrations and immersion times of the devices in the cleaning or disinfecting solutions indicated by the manufacturer.

2 Selection of material to carry out cleaning. For manual dirt removal, use only nylon brushes; Never use metal brushes or steel wool.

3 How to carry out cleaning/disinfection. Carefully remove all postoperative residues (Blood, bone...), in the case of products with cavities, such as drills with irrigation, rinse all cavities of the products using a disposable syringe. Scrub the devices with a soft-bristled nylon brush in tap water until all dirt has been removed. Immerse the devices in an ultrasonic bath with a cleaning solution at the concentration and temperature specified in the detergent manufacturer's instructions.

4 Rinse for one minute, preferably with purified or distilled water until the cleaning solution has been removed. If there are still traces of dirt after cleaning, repeat the same procedure.

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING CE 0051
	2.2.2 IFU (EXPANDER/ EXPANSORES)	Rev.: 03
		Date: 05/12/2024
		Página 11 de 13

Then, immerse in a bath with disinfectant solution at the concentration and temperatures indicated by the manufacturer. Make sure all ducts are filled with disinfectant solution and all surfaces are covered with the solution.

5 Rinse with purified or distilled water until all the disinfectant solution has been removed.

6 Dry the devices immediately after the cleaning and disinfection process with single-use wipes. Never store instruments while they are still wet.

D) How to perform sterilization

Moist heat sterilization should be performed in an autoclave at 134° (273°F) for 5 minutes. Make sure that the elements inside the autoclave are not rusty. Surgical material must be sterilized inside sterilization bags. Do not sterilize materials of different types in contact. The exposure temperature of the products and sterilization trays will not exceed 137°C.

Only for the United States: The validated parameters are Tª 270°F / 132°C for 4 minutes and a drying time of 20 minutes.

SURGICAL PROCEDURES

A) Instructions for use. Socket due to bone compaction (Figure 1):

Step 1: Drill a pilot socket using the 2mm pilot drill (not included in the set).

Step 2: Insert expander No. 1 until the threaded area disappears (about 8mm), rotating it slowly clockwise (right). Then remove it slowly by turning it in the opposite direction. **(*1)**

Step 3: Introduce expander No. 2 to the required depth depending on the implant to be inserted **(*2)** by rotating it slowly clockwise (clockwise). And then extract it slowly by turning it in the opposite direction. For 3.3 and 3.4mm diameter implants, skip to step 7.

Step 4: Continue the expansion sequence by slowly introducing expander No. 3 to the required depth. For 3.75mm diameter implants, skip to step 7.

Step 5: Slowly introduce expander No. 4 to the required depth. For 4.25mm diameter implants, skip to step 7.

Step 6: Complete the expansion sequence with expanders #5 and #6 to the depth required by the implant. **(*3)**

Step 7: Slowly insert the implant respecting the recommended insertion torque limit (35 - 45 Ncm). If necessary, stop for a few seconds before continuing with the insertion, or even go back if the insertion torque is excessive to prolong the expansion sequence until an appropriate insertion torque is achieved.

(*1) If excessive resistance is encountered during the introduction of any of the expanders, it is recommended to stop the expansion for a few seconds so as not to cause excessive heating before continuing to the desired depth.

(*2) See “Expander marking” and “Expander template” for selection of the required depth.

(*3) With the No. 6 dilator expander it is possible to perform sinus elevations via alveolar due to its convex tip.

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING CE 0051
	2.2.2 IFU (EXPANDER/ EXPANSORES)	Rev.: 03
		Date: 05/12/2024
		Página 12 de 13

Figure 1: Steps 1, 2, 3, 4 and 7 of the sequence of using the expanders to create the socket through bone compaction.

B) Instructions for use. Cortical Split with expanders (Figure 2):

Step 1b: Uncover the bone of the part to be treated, and remove the soft tissues, trying to maintain the vascularization of the area.

Step 2b: Make a longitudinal incision in the area to be treated using a surgical saw.

Step 3b: Mark, if necessary, the insertion point or points of the expanders using the 2mm pilot drill (not included in the set), or using a ball drill (not included in this set).

Step 4b: Follow the sequence described in steps 2 to 7, alternately inserting the corresponding expander into each of the marked expansion points.

Step 5b: Insert the implants, place the corresponding closure or healing screws and suture.


Figure 2: Main steps in the sequence of using the expanders for the “Cortical Split” performance.

PACKAGING AND STORAGE

The product is delivered clean, but not sterile. All information about the product is on the label. Product packaging is done individually. No special storage conditions are needed since the raw materials used are stable under normal pressure and temperature conditions. It is recommended to store between temperatures between 18°C and 50°C and a humidity limit below 50%. Store the set at a moderate temperature in a dry, dust-free and ventilated place where corrosive vapors are not produced.

DISPOSAL

The product should be disposed of in an environmentally friendly manner in accordance with local regulations. Hazardous waste of contaminated devices or sharp objects should be disposed of in suitable containers which meet the specific technical requirements.

DETECTION AND NOTIFICATION OF INCIDENTS WITH HEALTH PRODUCTS

Important notice: Any patient or healthcare personnel who detects any serious incident related to the Ticare dental implant must notify the manufacturer, Mozo Grau S.A. and the competent authority of the Member State in which the user and/or patient are established.

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING  0051
	2.2.2 IFU (EXPANDER/ EXPANSORES)	Rev.: 03
		Date: 05/12/2024
		Página 13 de 13

LABELLING SYMBOLS

SYMBOLS	DESCRIPTION
	Reference of the producto
	Batch code
	Manufacture date
	Manufacturer
	Medical device
	Non sterile
	CE Mark
	Consult electronic instructions for use
	Unique device identifier

MANUFACTURE

Mozo Grau S.A.
 C/ Santiago López González, N°7
 47197 Valladolid (Spain)
 Tel: +34 983 309 602
 e-mail: info@ticareimplants.com
www.ticareimplants.com