MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING  0051
	2.1.1 IFU (ABUTMENTS AND SCREWS CLASS IIb/ PROTESIS DENTALES CLASE IIb)	Rev.: 10
		Date: 17/12/2024
		Página 1 de 21

ESPAÑOL

ALCANCE

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ENGLISH

SCOPE

This instruction has the follow scope:

- Performance; Support piece for dental prosthesis.
 - o **Products: Straight abutments, Angled abutments, Tapered abutments, Titanium tapered abutment, Titanium Ucla abutments, Ball abutments.**
- Performance; piece with thread that will be used as joining element between the abutment and the implant or covering over the implant.
 - o **Products; Cover screws, Retaining screws, Healing abutments, Caps for tapered abutments**

This classification has been made based on performance.
These products are compatible with Inhex and Osseous Dental Implants.

DESCRIPTION

Prosthetic system is made of Grade V Titanium a biocompatible material. The prostheses are manufactured with High-precision machinery with a margin of error of ± 0.005 mm and a tolerance of ± 0.01 mm.

INTENDED USE

Medical device is intended for oral endosseous implantation in the upper and lower jawbone to provide support for dental prosthetic restorations of edentulous and partially dental patients.

INTENDER USER

TICARE Dental Implant System is intended to be used by qualified dental surgeons or dental healthcare professionals who have received additional training in dental implant procedures (dental implantology).

INTENDED TARGET GROUPS

TICARE Dental Implant System is to be used in edentulous and partially dentate patients, with stablished bone and finished growth period, and who are subject to dental implant treatments or aesthetic oral rehabilitation.

INDICATIONS (INCLUDING MEDICAL CONDITIONS)

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING  0051
	2.1.1 IFU (ABUTMENTS AND SCREWS CLASS IIb/ PROTESIS DENTALES CLASE IIb)	Rev.: 10
		Date: 17/12/2024
		Página 12 de 21

TICARE Dental Implant System is intended to be used by qualified dental surgeons of dental healthcare professionals who have received additional training in dental implant procedures. It is intended to be used in edentulous and partially dental patients, with established bone and who have finished its growth period.

The system is indicated for the following medical conditions:

- Oral endosseous implantation in the upper and lower jaw to restore chewing function.
- Aesthetic dental rehabilitation.

TICARE implants are indicated for single-stage or two-stage procedures to support screw retained restorations and can be used for immediate loading when good primary stability is achieved and with appropriate occlusal loading. Small diameter (3.3mm) TICARE Mini implants are intended to replace a lateral incisor in the maxilla and/or a central or lateral incisor in the mandible and should not be used in the molar region. TICARE Osseous Quattro and Inhex Quattro dental implants are indicated to support permanently fixed restorations and, due to their thread design, are not indicate for use in bone type I.

TICARE Inhex and Osseous implant of 6mm length are intended for use in a two-stage surgical procedure and are indicated for delayed loading to support permanently fixed restorations. These implants are indicated only for straight abutments.

Mozo Grau S.A. developed for each of their implant systems, various preprosthetic and prosthetic solutions to meet the particular needs the patients.

1) Cover Screws

Device Description: These abutments are for cover the Osseous implants during the osseointegration. To be placed on Osseous Implants. Healing abutments are available in different heights and diameters.

Material: Cover Screws are made of titanium grade 5, ASTM F136

Method of application: place immediately after implant placement using the 1.25 mm screwdriver. Apply manual torque (10 Ncm maximum).

Precautions: The platform of the implant should be clean from tissue or bone remnants which may prevent the proper seating of the abutment into the implant. Before placing the abutments, you should clean the inner thread of the implant.

2) Retaining Screws

Device Description: These screws are used for screwing the abutments to the implants and the restorations into the abutments.

Material: Prosthesis retention screws are made on titanium grade 5 according to ASTM 136.

Method of application: screw into the abutment or implant with the 1.25 mm hex screwdriver with the recommended torque.

Precautions: recommended tightening torques:

- Osseous direct to implant – 30 Ncm
- Inhex MINI direct to implant – 20 Ncm
- Inhex STD-MAXI direct to implant – 30 Ncm
- Osseous into tapered abutments – 20 Ncm
- Inhex MINI into tapered abutments – 15 Ncm
- Inhex STS-MAXI tapered abutments – 20 Ncm

3) Healing Abutments

A) Osseous Healing Abutments

Device Description: These abutments are used for modeling soft tissue. To be placed on Osseous Implants. Healing abutments are available in different heights and diameters.

Material: Healing abutments are made of titanium grade 5, ASTM F136

Method of application: place immediately after implant placement or in a delayed second surgery using the 1.25 mm screwdriver. Apply manual torque (10 Ncm maximum). It is convenient to leave the healing screw 2 mm protruding from the soft tissue.

For common healing situations of Mini, Standard and Maxi Osseous implants:

- Common situations for Mini Osseous implants: situations in which the space between two teeth is less than 4.5 mm.
- Common situations for Standard Osseous implants: situations in which the space between two teeth is between 4.5 and 5.5 mm.
- Common situations for Maxi Osseous implants: situations in which the space between two teeth is more than 5.5 mm.

Precautions: The platform of the implant should be clean from tissue or bone remnants which may prevent the proper seating of the abutment into the implant. Before placing the abutments, you should clean the inner thread of the implant.

B) Aesthetic Osseous Healing Abutments

Device Description: These abutments are used for modeling the soft tissue. To be placed on Osseous Implants. Healing abutments are available in different heights and diameters.

Material: Healing abutments are made of titanium grade 5, ASTM F136

Method of application: place immediately after implant placement or in a delayed second surgery using the 1.25 mm screwdriver. Apply manual torque (10 Ncm maximum). It is convenient to leave the healing screw 2 mm protruding from the soft tissue. For uncommon healing situations of Standard and Maxi Osseous implants:

- Uncommon situations for Standard Osseous implants: situations in which the space between two teeth is less than 4.5 and more than 5.5 mm.
- Uncommon situations for Maxi Osseous implants: situations in which the space between two teeth is less than 5.5 mm.

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING  0051
	2.1.1 IFU (ABUTMENTS AND SCREWS CLASS IIb/ PROTESIS DENTALES CLASE IIb)	Rev.: 10
		Date: 17/12/2024
		Página 13 de 21

Precautions: The platform of the implant should be clean from tissue or bone remnants which may prevent the proper seating of the abutment into the implant. Before placing the abutment, you should clean the inner thread of the implant.

C) InHex Healing Abutments

Device Description: These abutments are used for modeling the soft tissue. To be placed on InHex Implants. Healing abutments are available in different heights and diameters.

Material: Healing abutments are made of titanium grade 5, ASTM F136

Method of application: place immediately after implant placement or in a delayed second surgery using the 1.25 mm screwdriver. Apply manual torque (10 Ncm maximum). It is convenient to leave the healing screw 2 mm protruding from the soft tissue.

Precautions: The platform of the implant should be clean from tissue or bone remnants which may prevent the proper seating of the abutment on the implant. Before inserting the abutment, you should clean the inner thread of the implant.

D) 5mm InHex Healing Abutments

Device Description: These abutments are used for modeling soft tissue. To be placed on InHex Implants. Healing abutments are available in different heights and diameters.

Material: Healing abutments are made of titanium grade 5, ASTM F136

Method of application: place immediately after implant placement or in a delayed second surgery using the 1.25 mm screwdriver. Apply manual torque (10 Ncm maximum). For achieving a wider emergence profile. It is convenient to leave the healing screw 2 mm protruding from the soft tissue.

Precautions: The platform of the implant should be clean from tissue or bone remnants which may prevent the proper seating of the abutment on the implant. Before inserting the abutment, you should clean the inner thread of the implant.

E) Aesthetic InHex Healing Abutments

Device Description: These abutments are used for modeling soft tissue. To be place on InHex Implants. Healing abutments are available in different heights and diameters.

Material: Healing abutments are made of titanium grade 5, ASTM F136

Method of application: place immediately after implant placement or in a delayed second surgery using the 1.25 mm screwdriver. Apply manual torque (10 Ncm maximum). For achieving a wider emergence profile. It is convenient to leave the healing screw 2 mm protruding from the soft tissue.

Precautions: The platform of the implant should be clean from tissue or bone remnants which may prevent the proper seating of the abutment on the implant. Before inserting the abutment, you should clean the inner thread of the implant.

4) Healing and Provisional Abutments for Tapered Abutments

A) Healing Abutments for tapered abutments (one part)

Device Description: These abutments are used to cover the tapered abutments during the osseointegration. It is an abutment of a single part.

Material: Healing abutments for tapered abutments are made of titanium grade 5, ASTM F136

Method of application: Place on tapered abutments using the 1.25 mm hex screwdriver applying manual torque (maximum 10 Ncm)

Precautions: Be careful not pinching the soft tissue when placing it

B) Provisionals for tapered abutments (two parts)

Device Description: These abutments are used to cover the tapered abutments during the osseointegration or covering the tapered abutments in temporary restorations. They are composed of two parts: an abutment and a screw.

Material: Provisional abutments for tapered abutments are made of titanium grade 5, ASTM F136

Method of application: Place on tapered abutments using the 1.25 mm hex screwdriver applying manual torque (maximum 10 Ncm)

Precautions: Be careful not pinching the soft tissue when placing it. Make sure to release enough the existing denture in order to avoid overpressure on the implants.

5) Titanium Straight Abutments

A) Osseous Titanium Prepable Straight Abutments

Device Description: These abutments are designed for cement-retained single (hexagonal) or multiple-unit restorations (non hexagonal) direct to an Osseous implant. These abutments are intended for temporary restorations (for temporary abutments the maximum stay in the mouth is 2-3 months. The anatomical part of the tooth must be out of occlusion with the antagonist teeth to avoid transmitting loads to the temporary restoration). Soft tissue height must be above the abutment transmucosal ring.

Material: Titanium Straight abutments are made of titanium grade 5, ASTM F136

Method of application: Place the abutment already prepared for the lab technician into the implant and screw it at 30 Ncm using the 1.25 mm hex screwdriver. Impression taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

Warning: Only modification available to the user is to shorten the post height to a minimum of 4mm. No other modification is recommended.

B) Osseous Titanium Straight Abutments

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING  0051
	2.1.1 IFU (ABUTMENTS AND SCREWS CLASS IIb/ PROTESIS DENTALES CLASE IIb)	Rev.: 10
		Date: 17/12/2024
		Página 14 de 21

Device Description: These abutments are designed for cement-retained single (hexagonal) or multiple-unit restorations (non hexagonal) direct to an Osseous implant. These abutments are designed for the general cases of definitive restorations. Soft tissue height must be above the abutment transmucosal ring.

Material: Titanium Straight abutments are made of titanium grade 5, ASTM F136

Method of application: Place the abutment already prepared for the lab technician (it already has a shoulder but needs to be adapted to the specific case) into the implant and screw it at 30 Ncm using the 1.25 mm hex screwdriver. Impression taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

C) Osseous Titanium Aesthetic Straight Abutments

Device Description: These abutments are designed for definitive cement-retained single- or multiple-unit restorations on Osseous implants. These abutments are designed for cases in which the restorative space is 6 mm or more or a wide emergence profile is desired.

Material: Titanium Straight abutments are made of titanium grade 5, ASTM F136

Method of application: Place the abutment already prepared for the lab technician (it already has a shoulder but needs to be adapted to the specific case) into the implant and screw it at 30 Ncm using the 1.25 mm hex screwdriver.

Impression taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

D) Osseous Titanium Multitask Abutments

Device Description: These abutments are designed for definitive cement-retained single- or multiple-unit restorations on Osseous implants.

Material: Titanium Straight abutments are made of titanium grade 5, ASTM F136

Method of application: Place the abutment already prepared for the lab technician (it already has a shoulder but needs to be adapted to the specific case) into the implant and screw it at 30 Ncm using the 1.25 mm hex screwdriver. Impression taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

Warning: Only modification available to the user is to shorten the post height to a minimum of 4mm. No other modification is recommended.

E) InHex Titanium Preapable Straight Abutments

Device Description: These abutments are designed for cement-retained single (hexagonal) or multiple-unit restorations (no hexagonal) direct to an InHex implant. These abutments are intended for temporary restorations (for temporary abutments the maximum stay in the mouth is 2-3 months The anatomical part of the tooth must be out of occlusion with the antagonist teeth to avoid transmitting loads to the temporary restoration). Soft tissue height must be above the abutment transmucosal ring.

Material: Inhex Titanium abutments are made of titanium grade 5, ASTM F136

Method of application: Place the abutment already prepared for the lab technician into the implant and screw it at 30 Ncm using the 1.25 mm hex screwdriver. Impression taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

Warning: Only modification available to the user is to shorten the post height to a minimum of 4mm. No other modification is recommended.

F) InHex Titanium Straight Abutments

Device Description: These abutments are designed for definitive cement-retained single- or multiple-unit restorations on Inhex implants.

Material: InHex Titanium Straight abutments are made of titanium grade 5, ASTM F136

Method of application: Place the abutment already prepared for the lab technician (it already has a shoulder but needs to be adapted to the specific case) into the implant and screw it at 20 (MINI) or 30 (STD and MAXI) Ncm using the 1.25 mm hex screwdriver. Impression taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

G) InHex Titanium Aesthetic Straight Abutments

Device Description: These abutments are designed for definitive cement-retained single- or multiple-unit restorations on InHex implants. These abutments are designed for cases in which the restorative space is 6 mm or more or a wide emergence profile is desired.

Material: InHex Titanium Straight abutments are made of titanium grade 5, ASTM F136

Method of application: Place the abutment already prepared for the lab technician (it already has a shoulder but needs to be adapted to the specific case) into the implant and screw it at 30 Ncm using the 1.25 mm hex screwdriver.

Impression taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

H) InHex Titanium 5mm Straight Abutments

Device Description: These abutments are designed for definitive cement-retained single- or multiple-unit restorations on InHex implants. These abutments are designed for cases in which the restorative space is 5 to 6 mm or more or a wide emergence profile is desired.

Material: InHex Titanium Straight abutments are made of titanium grade 5, ASTM F136

Method of application: Place the abutment already prepared for the lab technician (it already has a shoulder but needs to be adapted to the specific case) into the implant and screw it at 30 Ncm using the 1.25 mm hex screwdriver.

Impression taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

I) Inhex Titanium Multitask Abutments

Device Description: These abutments are designed for definitive cement-retained single- or multiple-unit restorations on Inhex implants.

Material: InHex Titanium Straight abutments are made of titanium grade 5, ASTM F136

Method of application: Place the abutment already prepared for the lab technician (it already has a shoulder but needs to be adapted to the specific case) into the implant and screw it at 20 (MINI) or 30 (STD and MAXI) Ncm using the 1.25 mm hex screwdriver. Impression

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING  0051
	2.1.1 IFU (ABUTMENTS AND SCREWS CLASS IIb/ PROTESIS DENTALES CLASE IIb)	Rev.: 10
		Date: 17/12/2024
		Página 15 de 21

taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

Warning: Only modification available to the user is to shorten the post height to a minimum of 4mm. No other modification is recommended.

6) Titanium Angled Abutments

A) Osseous Titanium Angled Abutments

Device Description: These abutments are designed for definitive cement-retained single or multiple-unit restorations into an Osseous implant. These abutments are designed for cases in which the emergence angulation of the crown or disparallelism between implants must be corrected.

Material: Osseous Titanium angled abutments are made of titanium grade 5, ASTM F136.

Method of application: Place the abutment already prepared for the lab technician (it already has a shoulder but needs to be adapted to the specific case) into the implant and screw it at 30 Ncm using the 1.25 mm hex screwdriver. Impression taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

B) InHex Titanium Angled Abutments

Device Description: These abutments are designed for cement-retained single or multiple-unit restorations direct to an InHex implant. These abutments are designed for the general cases of definitive restorations in which the emergence angulation of the crown or disparallelism between implants must be corrected.

Material: InHex Titanium angled abutments are made of titanium grade 5, ASTM F136.

Method of application: Place the abutment already prepared for the lab technician (it already has a shoulder but needs to be adapted to the specific case) into the implant and screw it at 20 (MINI) or 30 (STD and MAXI) Ncm using the 1.25 mm hex screwdriver. Impression taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

C) InHex Titanium 5mm Angled Abutments

Device Description: These abutments are designed for cement-retained single or multiple-unit restorations direct to an InHex implant. These abutments are designed for the general cases of definitive restorations in which the emergence angulation of the crown or disparallelism between implants must be corrected and in addition the restorative space is more than 5 mm or a wider emergence profile is desired.

Material: InHex Titanium angled abutments are made of titanium grade 5, ASTM F136.

Method of application: Place the abutment already prepared for the lab technician (it already has a shoulder but needs to be adapted to the specific case) into the implant and screw it at 20 (MINI) or 30 (STD and MAXI) Ncm using the 1.25 mm hex screwdriver. Impression taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

7) Titanium Tapered Abutments

A) Osseous Titanium Tapered Abutments

Device Description: These abutments are designed for screw-retained multiple-unit restorations on Osseous Implants and for overdentures. They allow to correct disparallelism between implants up to 60°.

Material: Osseous Titanium Tapered abutments are made of titanium grade 5, ASTM F136.

Method of application: Place the abutment directly into the implants and tighten it at 30 Ncm. Impression taking at abutment-level. Lab work with abutment analogs. Screw the restoration at 20 Ncm.

Precautions: Recommended for soft tissue height above 3 mm

B) Osseous Titanium Straight Tapered Abutments

Device Description: These abutments are designed for screw-retained multiple-unit restorations on Osseous Implants and for overdentures. They allow to correct disparallelism between implants up to 40°.

Material: Osseous Titanium Tapered abutments are made of titanium grade 5, ASTM F136.

Method of application: Place the abutment directly into the implants and tighten it at 30 Ncm. Impression taking at abutment-level. Lab work with abutment analogs. Screw the restoration at 20 Ncm.

Precautions: Recommended for soft tissue height above 3 mm.

C) Osseous Titanium Angled Tapered Abutments

Device Description: These abutments are designed for screw-retained multiple-unit restorations on Osseous Implants and for overdentures. They allow to correct disparallelism between implants up to 60°. They allow to correct the positioning of the screw hole into a more suitable place.

Material: Osseous Titanium Tapered abutments are made of titanium grade 5, ASTM F136.

Method of application: Place the abutment directly into the implants and tighten it at 30 Ncm. Impression taking at abutment-level. Lab work with abutment analogs. Screw the restoration at 20 Ncm.

Precautions: Recommended for soft tissue height above 3 mm.

D) InHex Titanium Tapered Abutments

Device Description: These abutments are designed for screw-retained multiple-unit restorations on Inhex Implants and for overdentures. They allow to correct disparallelism between implants up to 60°.

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING  0051
	2.1.1 IFU (ABUTMENTS AND SCREWS CLASS IIb/ PROTESIS DENTALES CLASE IIb)	Rev.: 10
		Date: 17/12/2024
		Página 16 de 21

Material: Inhex Titanium Tapered abutments are made of titanium grade 5, ASTM F136.

Method of application: Place the abutment directly into the implants and tighten it at 20 (MINI) or 30 (STD and MAXI) Ncm. Impression taking at abutment-level. Lab work with abutment analogs. Screw the restoration at 15 (MINI) or 20 (STD and MAXI) Ncm.

Precautions: Recommended for soft tissue height above 3 mm

E) InHex Titanium Straight Tapered Abutments

Device Description: These abutments are designed for screw-retained multiple-unit restorations on Inhex Implants and for overdentures. They allow to correct disparallelism between implants up to 40°.

Material: Inhex Titanium Tapered abutments are made of titanium grade 5, ASTM F136.

Method of application: Place the abutment directly into the implants and tighten it at 20 (MINI) or 30 (STD and MAXI) Ncm. Impression taking at abutment-level. Lab work with abutment analogs. Screw the restoration at 15 (MINI) or 20 (STD and MAXI) Ncm.

Precautions: Recommended for soft tissue height above 3 mm

F) InHex Titanium Angled Tapered Abutments

Device Description: These abutments are designed for screw-retained multiple-unit restorations on Inhex Implants and for overdentures. They allow to correct disparallelism between implants up to 60°. They allow to correct the positioning of the screw hole into a more suitable place.

Material: Inhex Titanium Tapered abutments are made of titanium grade 5, ASTM F136.

Method of application: Place the abutment directly into the implants and tighten it at 30 Ncm. Impression taking at abutment-level. Lab work with abutment analogs. Screw the restoration at 20 Ncm.

Precautions: Recommended for soft tissue height above 3 mm

8) Titanium Ball Abutments

Osseous/Inhex Titanium Ball Abutments

Device Description: These abutments are designed for attachment-retained restoration on Osseous Implants, 10° of maximum angulation between implants.

Material: Titanium ball abutments are made of titanium grade 5, ASTM F136.

Method of application: Impression taking at abutment level. Place it into the implant at 30 Ncm.

Precautions: It is convenient that the transmucosal part of the abutment emerges at least 1 mm from the gingiva.

9) UCLA Abutments

A) Osseous UCLA Abutments

Device Description: For cement-retained restorations, single or multiple, direct on Osseous implants in challenging situations. Also, for welding titanium bars.

Material: Osseous UCLA abutments are made of titanium grade 5, ASTM F136.

Method of application: Place the abutment already prepared for the lab technician into the implant and screw it at 30 Ncm using the 1.25 mm hex screwdriver. Impression taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

B) Osseous Temporary UCLA Abutments

Device Description: For screw-retained temporary restorations, single or multiple, direct on Osseous implants. The anatomical part of the tooth must be out of occlusion with the antagonist teeth to avoid transmitting loads to the temporary restoration. The maximum stay in the mouth is 2-3 months.

Material: Osseous UCLA abutments are made of titanium grade 5, ASTM F136.

Method of application: abutment that allows chair- or lab-side preparation. Screw it into the implant at 10 Ncm using the 1.25 mm hex screwdriver.

Precautions: Temporary abutment.

C) InHex UCLA Abutments

Device Description: For cement-retained restorations, single or multiple, direct on Osseous implants in challenging situations. Also, for welding titanium bars.

Material: Inhex UCLA abutments are made of titanium grade 5, ASTM F136.

Method of application: Place the abutment already prepared for the lab technician (it already has a shoulder but needs to be adapted to the specific case) into the implant and screw it at 30 Ncm using the 1.25 mm hex screwdriver. Impression taking on implant level.

Precautions: Soft tissue height should be slightly above the abutment shoulder.

D) InHex Temporary UCLA Abutments

Device Description: For screw-retained single-unit temporary restorations, direct on Inhex implants. The anatomical part of the tooth must be out of occlusion with the antagonist teeth to avoid transmitting loads to the temporary restoration. The maximum stay in the mouth is 2-3 months.

Material: Inhex UCLA abutments are made of titanium grade 5, ASTM F136.

Method of application: abutment that allows chair- or lab-side preparation. Screw it into the implant at 10 Ncm using the 1.25 mm hex screwdriver.

Precautions: Temporary abutment.

10) Titanium Posts for tapered Abutments

A) Temporary Titanium Pots for Tapered Abutments

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING CE 0051
	2.1.1 IFU (ABUTMENTS AND SCREWS CLASS IIb/ PROTESIS DENTALES CLASE IIb)	Rev.: 10
		Date: 17/12/2024
		Página 17 de 21

Device Description: These abutments are designed for temporary screw-retained multiple-unit restorations on Osseous or InHex tapered abutments. The anatomical part of the tooth must be out of occlusion with the antagonist teeth to avoid transmitting loads to the temporary restoration. The maximum stay in the mouth is 2-3 months.
 Material: Titanium posts for tapered abutments are made on titanium grade 5 according to ASTM F136.
 Method of application: preparation at chair- or lab-side. Place it on tapered abutments with the 1.25 hex screwdriver at 10Ncm.
 Precautions: For temporary restorations


B) Titanium Pots for Tapered Abutments
 Device Description: These posts are designed for screw-retained multiple-unit restorations on Osseous or InHex tapered abutments in challenging situations. Also, for welding titanium bars.
 Material: Titanium posts for tapered abutments are made on titanium grade 5 according to ASTM F136.
 Method of application: Screw the post into the abutment at 15 (MINI) or 20 (STD and MAXI) NCM using the 1.25 mm hex screwdriver.
 Precautions: Soft tissue height should be slightly above the abutment shoulder.

WARNING; Abutments with dimensional constraints

1) Osseous Titanium Straight Abutments

Min: It is 4 mm for every Osseous Prepable Straight Abutment, Osseous Straight Abutment, Osseous Aesthetic Straight Abutment and Osseous Multitask Abutment.


Max: It is 7 mm for every Osseous Prepable Straight Abutment, Osseous Straight Abutment, Osseous Aesthetic Straight Abutment and 5.10 mm for every Osseous Multitask Abutment.


2) InHex Titanium Straight Abutments

Min: It is 4 mm for every Inhex Prepable Straight Abutment, InHex Straight Abutment, InHex Aesthetic Straight Abutment, Inhex 5 mm Straight Abutment and InHex Multitask Abutment.

Max: It is 6.5 mm for every Inhex Prepable Straight Abutment, InHex Straight Abutment, InHex Aesthetic Straight Abutment, Inhex 5 mm Straight Abutment and 5.10 mm for InHex Multitask Abutment.

3) Osseous UCLA Abutments


For every Osseous UCLA:

R1(mm)	R2 (mm)	R3 (mm)	R5
4	9	0.80	0°

(R4 and R6 are determined by the morphology of the mouth in the rehabilitation area)

4) InHex UCLA Abutments


For every InHex UCLA:

R1(mm)	R2 (mm)	R3 (mm)	R5
4	9	Between 0.60 y 0.70	0°

(R4 and R6 are determined by the morphology of the mouth in the rehabilitation area)

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING CE 0051
	2.1.1 IFU (ABUTMENTS AND SCREWS CLASS IIb/ PROTESIS DENTALES CLASE IIb)	Rev.: 10
		Date: 17/12/2024
		Página 18 de 21

5) Titanium Posts for Tapered Abutments

For every Post for 10° Tapered Abutment and Post for Angled Tapered Abutment

R1(mm)	R2 (mm)	R3 (mm)	R5
5	7.9	0.80	0°

(R4 and R6 are determined by the morphology of the mouth in the rehabilitation area)

For every Post for 30° Tapered Abutment

R1(mm)	R2 (mm)	R3 (mm)	R5
4	7.9	0.80	0°

(R4 and R6 are determined by the morphology of the mouth in the rehabilitation area)

CONTRAINDICATIONS

Lack of quality and quantity of remaining bone, infection, and local or systemic factors that may interfere with the healing process of bone or soft tissue (poor oral hygiene, uncontrolled diabetes, pregnancy, oral radiotherapy, i.a.) are some potential causes of stability loss.

Dental implant procedures are contraindicated for patients who are medically unfit to undergo an oral surgical procedure.

- Patients allergic to titanium alloy materials.
- Patients with insufficient mental health precluding patient cooperation.
- Patients with abuse of drugs and alcohol
- Patients who have suffered myocardial infarct last year, oral infections, or malignancies
- Patients with uncontrolled diabetes or blood disorders

Information to allergen

Do not use if the patient presents known allergy or hypersensitivity to any of the material components.

- Abutments; Titanium Gr. 5 (Ti, Al, V, Fe, O, N, C)

GENERAL WARNINGS

- Do not use if the individual package is open or damaged.
- To be used by qualified professionals only (dentists, doctor experts in dentistry or stomatology, maxillofacial surgeons).
- For the safety and effectiveness of use of Inhex / Osseous Prosthetic System, it is strongly suggested that specific training be carried out, as well as the study of appropriate textbooks. THESE INSTRUCTIONS ARE NOT INTENDED TO BE SUBSTITUTE FOR ADEQUATE TRAINING. Responsibility for the proper choice of the patient, the adequate training, the experience in implant placement of the implants, and provision of appropriate information for informed consent devolves upon the practitioner. Improper technique may result in implant failure and/or loss of supportive bone. Mozo Grau is not able to control the application and processing of our products, which are user's responsibility. Therefore, Mozo Grau will not accept liability for damage arising thereof.

SPECIFIC WARNINGS

- In Treatment Planning Prior to Surgery:
 - o Warning: It is highly recommendable to clean, disinfect and sterilize all the Inhex / Osseous Prosthetic tools before each use, specially before first use after delivery since products are supplied non-sterile (they must be cleaned and disinfected after having removed packaging).
A proper cleaning and disinfection are essential prerequisites for proper sterilization.
Please do not forget that the disinfectant used in the pre-treatment is only for personal protection and does not replace the later disinfection after cleaning. See section CLEANING AND STERILIZATION.
 - o Warning: mechanical fracture of implants or bridges may occur. A blow to the jaw or mouth, or the stress concentration coming from the bridge, may result in such mechanical failures. Do not use small diameter implants in high load applications. Mobility, bone loss, and / or pain and infection may indicate that the implant is failing. On rare occasions when the implant

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING  0051
	2.1.1 IFU (ABUTMENTS AND SCREWS CLASS IIb/ PROTESIS DENTALES CLASE IIb)	Rev.: 10
		Date: 17/12/2024
		Página 19 de 21

osseointegration fails, the implant has to be removed, then, the bone heals and the condition of the jaw and mouth is virtually the same as if the procedure had not been performed. The loss of implant stability (osseointegration failure) and loss of the bridge are possible incidents after surgery.

- Careful clinical and radiological examination of the patient is essential to treatment planning prior to surgery.

- In Handling:

- Warning: No Sterile, single-use device. Do not re-use. Do not re-process. Re-use could cause cross contamination. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Use of a non-sterile device may lead to infection. Titanium is a metal which may be easily scratched during use. A scratched surface is more difficult to clean existing, therefore, a greater risk of contamination.
- Warning: Recommended torque level is 30Ncm for abutment attaches directly to implant, 20Ncm for transepithelial abutments and healing screws 10Ncm

GENERAL PRECAUTIONS

- Caution: Federal (US) law restricts this device to sale by or on the order of a licensed dentist or physician.
- Ticare dental abutments are intended to be used only in conjunction with Ticare dental implants and associated instruments according to the instructions and recommendations of Mozo Grau. Combining implants, instruments and prosthetic components that are not dimensioned for correct mating can lead to mechanical failure, damage to tissue or unsatisfactory aesthetic results.

SPECIFIC PRECAUTIONS

- In Treatment Planning Prior to Surgery
 - Precaution: Special care is advised when treating patients who have local or systemic risk factors that may interfere with the healing process of bone or soft tissue. These risk factors include smoking, poor oral hygiene, uncontrolled diabetes, pregnancy, oral radiotherapy, infections in neighboring bone and certain therapies including steroids and chemotherapy. Bruxism or other jaw disorder may impede the success of dental implant procedures and may require the consideration of alternative treatment options.
 - Advise patients who smoke to stop smoking before undergoing dental implant procedures.
 - Dental implant procedures are associated with auxiliar medical acts (ex: x-rays, administration of anaesthetics, antibiotic treatments, etc.) that may be contra-indicated in some medical conditions such as pregnancy or breastfeeding amongst others. It is the final responsibility of the dental practitioner or qualified professional to judge whether the patients are medically fit or unfit to undergo an oral surgical procedure
 - Regular dental check-ups and maintaining good oral hygiene are crucial to minimize any adverse effects..
 - Routine use of abutments and screws for dental implants in paediatric patients is not recommended until jaw bone growth is completed.
- In Handling and Sterilization
 - Precaution: All instruments used should be sterilized before use. Use gloves to handle the product. Avoid dust and contamination of the surgical site. See section CLEANING AND STERILIZATION.
- In Directions for Use
 - PRECAUTION: Dental abutments are small in size. Care must be taken to ensure that abutments are not swallowed or aspirated by the patient.
 - PRECAUTION: Avoid overtightening the implant. Do not exceed the recommended maximum torque. Overtightening the abutments may to damage or fracture the screws.

CLEANING AND STERILIZATION

It is recommendable to clean, and sterilize all the Inhex / Osseous Prosthetic before use after delivery since products are supplied non-sterile (they must be cleaned and disinfected after having removed packaging).

A proper cleaning and disinfection are essential prerequisites for proper sterilization.

Recommended cleaning method

- Immerse devices in cleaning solution.
 - Mozo Grau does not recommended any specific cleaning and/or disinfection agents. Only specifically formulated cleaning agents should be used. Solution with a pH value of 7-10 and a temperature of approx. 35°C/95°C.
- Rinse and flush with purified or sterile water. Mineral residues from hard water can result in staining of the device.
- Dry with compressed air or wipes.

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING  0051
	2.1.1 IFU (ABUTMENTS AND SCREWS CLASS IIb/ PROTESIS DENTALES CLASE IIb)	Rev.: 10
		Date: 17/12/2024
		Página 20 de 21

Please do not forget that the disinfectant used in the pre-treatment is only for personal protection and does not replace the later disinfection after cleaning.

Method of sterilization

- Place product inside sterilization bags.
- Place the product inside the autoclave. Do not sterilize different types of materials in contact.
- Sterilization in autoclave by moisty heat 134 ° C (273 ° F) - 5 minutes. Make sure elements inside the autoclave are not rusty. Only for the United States: Validated parameters are Temperature 270°F / 132°C for 4 minutes and a drying time of 20 minutes.
- The exposure temperature for products and sterilization trays shall not exceed 137 ° C.
- When the autoclave is finished, remove the product from the autoclave.

SIDE EFFECTS AND POST OPERATIVE PROCEDURES IN THE SURGICAL AND PROSTHETIC PROCEDURES

For detailed information of the surgical and prosthetic procedures refer to the Surgical and Prosthetic manuals.

- In Side-Effects and Post-Operative Procedures
 - **PRECAUTION:** Avoid infection in the post-operative phase. Advise the patient to maintain strict oral hygiene. Smoking increases the risk of infection. Nicotine may interfere with the healing process. Advise the patient to stop smoking when undergoing dental implant procedures.
 - **PRECAUTION:** Regular patient follow-up is needed to secure long-term treatment outcome and to inform the patient about appropriate oral hygiene.

PACKAGING AND STORAGE

Prosthetic system is supplied clean but not sterile. All information about the product is labeled. Product packaging is performed individually.

This product does not need special storage and transport conditions because raw material is stable under pressure and atmospheric temperature conditions. It is recommended to store the product at temperatures under 50°C and humidity limit up to 70%

DISPOSAL

The product should be disposed of in an environmentally friendly manner in accordance with local regulations. Hazardous waste of contaminated devices or sharp objects should be disposed of in suitable containers which meet the specific technical requirements.

DETECTION AND NOTIFICATION OF INCIDENTS WITH HEALTH PRODUCTS

Important notice: Any patient or healthcare personnel who detects any serious incident related to the Ticare dental implant must notify the manufacturer, Mozo Grau S.A. and the competent authority of the Member State in which the user and/or patient are established.

LEGAL COMPLIANCE

After the launch of the new European Database on Medical Devices (Eudamed), the Summary of Safety and Clinical Performance (SSCP) will be available on the Eudamed website at <https://ec.europa.eu/tools/eudamed>. Until that time, you can find the SSCP at <https://www.ticareimplants.com/calidad/>

IMPLANT CARD

Mozo Grau S.A., as a manufacturer of implantable medical devices, is obliged to provide an **Implant card** with each dental implant, in accordance with Article 18 of the Medical Device Regulation (EU) 2017/745, which includes all the information necessary to identify both the implanted device and the patient and clinic where the implantation surgery has taken place.

LABELLING SYMBOLS

The following symbols are used on the packaging labels:

SYMBOLS	DESCRIPTION
	Reference of the producto
	Batch code

MOZOGRAU®	Technical File. Section 2: INFORMATION TO BE SUPPLIED BY THE MANUFACTURER	CE MARKING 
	2.1.1 IFU (ABUTMENTS AND SCREWS CLASS IIb/ PROTESIS DENTALES CLASE IIb)	Rev.: 10
		Date: 17/12/2024
		Página 21 de 21

	Manufacture date
	Manufacture
	Medical device
	Do not reuse
	Non esterile
	CE Mark
	Consult electronic instructions for use
	Unique device identifier

MANUFACTURE

Mozo Grau S.A.
C/ Santiago López González, N°7
47197 Valladolid (Spain)
Tel: +34 983 309 602
e-mail: info@ticareimplants.com
www.ticareimplants.com