ALCANCE

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MARCADO CE

SCOPE

This instructions has the follow scope: Fixing screw: Clase IIb

- 14101003 1.0 x 3.0mm Fixing screws
- 14101005 1.0 x 5.0mm Fixing screws
- 14101007 1.0 x 7.0mm Fixing screws
- 14101307 1.3 x 7.0mm Fixing screws
- 14101309 1.3 x 9.0mm Fixing screws
- 14101311 1.3 x 11.0mm Fixing screws
- 14101313 1.3 x 13.0mm Fixing screws

Fixing screws are used with the following Surgical devices class I:

- 14100563 1.0mm Fixing screw drill
- 14100569 1.3mm Fixing screw drill
- 14100564 Drill for boneblock
- 14100565 1.0mm Screwdriver tip
- 14100558 1.3mm Screwdriver tip
- 14100559 Screwdriver handpiece
- 14100566 Sterilization Box

*Titanium Mesh not included

INDICATIONS

A simple, fast and convenient way to fix bone grafts with total reability during the surgeries.

CONTRAINDICATIONS

Contraindications include:

1. Patients allergic to titanium alloy.
2. Patients with insufficient mental health that prevent patient cooperation
3. Patients who abuse drugs and alcohol

WARNINGS

For the security and effectiveness of the use of the system it is strongly suggested that a specific training expected to be undertaken, as well as the study of appropriate books. THESE INSTRUCTIONS ARE NOT INTENDED TO BE SUBSTITUTE FOR ADEQUATE TRAINING.

Information to allergen:

The materials used to manufacture the abutments are:

- Titanium Gr.5 (Ti, Al, V, Fe, O, N, C)

PRECAUTIONS

Fixing screws are for single-use. Re-using fixing screws may lead to infection and/or cross contamination. Fixing screws have to be used in aseptic conditions.

CLEANING AND STERILIZATION

Before use, clean and sterilise the Bone Grafting Set. The Bone Grafting Set should be cleaned, disinfected and sterilized before use. The products are supplied non-sterile. A proper cleaning is essential for proper sterilization.

Recommended cleaning method:

- Immerse devices in cleaning solution. Mozo Grau does not recommend any specific cleaning and/or disinfection agents. Only specifically formulated cleaning agents should be used. Solution with a pH value of 7-10 and a temperature of approx. 35°C/95°C.
- Rinse and flush with purified or sterile water. Mineral residues from hard water can result in staining of the device.
- Dry with compressed air or wipes.

Please do not forget that the disinfectant used in the pre-treatment is used for personal protection only, and it does not replace the disinfection to be carried out later, once the cleaning is performed.

Method of sterilization:

- Put the product inside the sterilization bags
- Place the product inside of the autoclave. Do not sterilize materials of different types in contact.
- Autoclave moist heat sterilization 134°C (237°F) – 5 minutes. Make sure that the elements inside the autoclave from rusted.
- The exposure temperature and product sterilization trays shall not exceed 137°C.
- When the autoclave is finished, remove the product from inside the autoclave.

SURGICAL PROCEDURES

One the donor area has been chosen, and before proceeding to the osteotomy, perform of holes must be done with the drill for bone blocks graft.

Next, remove the block graft by using chisels and then overlap it in the receiving area and through the holes previously made and in this case, with the drill for the screw of the chosen diameter. Perform the perforations in the receiving area.

Once this step is done, adjust the screws through the graft to the bone.

SEQUENCE OF USE OF BONE GRAFTING SET

Marking with 1.8mm
Lindemann Drill
(Ref. 13716718)

Cutting bone graft
perimeter with 1.8mm
Lindemann Drill
(Ref. 13716718)

Hole drilling in the graft
block with the Graft Block
Drill
(Ref. 14100564)

Releasing of Graft Block
with the Curved
Osteotome
(Ref. 14001031)

Receptor site drilling through the
Graft Block with the Drill for 1.0mm
screw (Ref. 14100563) or Drill for
1.3mm (Ref. 14100569) depending
on the size of the Graft Block.

Screwing the corresponding
screw (1.0 or 1.3mm)
different heights depending
on the case

Graft Block properly
fixed to the receptor site

Surgery site suture
covering both the donor
and the receptor site

PACKAGING AND STORAGE

Prosthetic system is supplied clean but not sterile. All information about the product is labeled. Product packaging is performed individually. Devices should be restored at room temperature.

This product does not need special conditions for storage and transport, because the raw material is stable under pressure and atmospheric conditions. It is recommended to store the product at a temperature between 10°C < T < 45°C (bags) or T < 50°C (Bilster).

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REFERENCE

LOT

MANUFACTURER

ONE USE ONLY, NO
REUSABLE

PAY ATTENTION TO
INSTRUCTIONS FOR USE

CE MARKING

0051