
1-Application

Ticare MultiUnit Straight Abutments are transgingival abutments that allow maintaining the gingival and bone level, improving the emergence profile and protecting the peri-implant tissues, facilitating the hygiene of the prosthesis. Ticare MultiUnit Straight Abutments allow to screw multiple restorations in gums higher than 3 mm. Due to the 20° of the cone angle, they help to correct disparelisms between implants, close to 40°.

2-Instructions for use


- Step 1:** Remove the set made of the MultiUnit Straight Abutment and the peek implant mount from its package and take it manually to the mouth, holding it by the peek piece.
- Step 2:** Manually screw it onto the implant until the abutment cone contacts the implant cone.
- Step 3):** Remove the peek implant mount from the abutment laterally displacing the implant mount softly (Image 1)

Step 4: Finally, use a ratchet wrench with a ratchet extender (or adapter to C / A direct A / C) to tighten the MultiUnit Straight Abutment onto the implant with a torque of 30 Ncm.

For the correct tightening of the MultiUnit Straight Abutment

- 1) Match the internal hexagon of the ratchet extender (or C / A adapter) with the external hexagon of the MultiUnit Straight Abutment cone (Image 2),
- 2) Match both part platforms and apply torque pressing slightly the set from the top of the ratchet wrench (Image 3).

Contraindications

For the safety and effectiveness of the instrumental this product is strongly suggested to be used by qualified personnel. THESE INSTRUCTIONS ARE NOT INTENDED TO BE A SUBSTITUTE FOR A PROPER TRAINING. This product is supplied NON-STERIL and must be used under aseptic conditions.

Warnings And Cautions

Those products whose characteristics have been altered by use and are not safe for the indicated functions shall not be used.

Product Composition

The Multi Unit abutments are made of titanium and the special wrench is made of peek

Packaging And Storage

All the product information is shown on the label. Product packaging is performed individually. No special storage conditions are needed as the raw materials used are stable under standard pressure and temperature conditions. It is advisable to store between 18°C and 50°C and a moisture limit below 70%. Store the ratchet at a well-ventilated place, in which there is no corrosive steam.

Cleaning And Sterilization

The products are supplied non-sterile. A proper cleaning and disinfection are prerequisites for proper sterilization.. Please do not forget that the disinfectant used in the pre-treatment used for personal protection only and does not replace the disinfection to be carried out later, one executed cleaning.

Method od sterilization:

- Put the producto inside retort pouches.
- Put the producto insided of the autoclave. Do not sterilize different types of materials in contact.
- Sterilization must be carried out in autoclave wet heat 134°C (237°F) – 5 minutes. Make sure that the elements inside the autoclave from corrosion.
- The exposure temperatura and product sterilization trays shall not exceed 137°C.
- When the autoclave is complete, remove the product inside the autoclave.

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