




INHEX PROSTHETIC SYSTEM

 ticare
i n H E X



GAP **ZERO**
T E C H N O L O G Y



Torque guide

		MINI	STD	MAXI
Cover screw		10 Ncm	10 Ncm	10 Ncm
Healing screw		10 Ncm	10 Ncm	10 Ncm
Impression coping		10 Ncm	10 Ncm	10 Ncm
Scanbody		3 Ncm	3 Ncm	3 Ncm
Customized abutments		20 Ncm	30 Ncm	30 Ncm
Final abutments (direct to the implant)		20 Ncm	30 Ncm	30 Ncm
Prosthesis over 10° or 30° tapered abutment		15 Ncm	20 Ncm	20 Ncm
Prostheses over 20° abutments (straight and angled multiunits)			15 Ncm	15 Ncm
Temporary restoration		10 Ncm	10 Ncm	10 Ncm
Screws for angled chimmney		Apply to the torque recommended by the manufacturer requierd for straight chimney		

Torque Ncm

MINI Platform restorations

HEALING ABUTMENTS



23203320 2 mm
23203330 3 mm
23203340 4 mm

IMPRESSION COPINGS



[Open tray]

23203303 Standard



23235504 Short L=14mm
23235502 Long L=20mm
23235505 Wax L=29



[Hex] [Closed tray]

23203324 Short
23203334 Long

[Non Hex] [Closed tray]
23205053

TCP

HEALING ABUTMENT



23603320 2 mm
23603330 3 mm

IMPRESSION COPINGS



[Open tray]

23603303



[Close tray]

23603324 Short
23603334 Long



SCANBODY

[Implant-level]

41636001 TC10

SCREW-RETAINED RESTORATIONS

SINGLE



BioCAM

TITANIUM INTERFACE

[Hex]

23237905 0,5 mm
23237910 1 mm
23237920 2 mm
23237930 3 mm



TCP

23637920 2 mm
23637930 3 mm

SCANBODY



[Implant-level]
41236001

CEMENT-RETAINED RESTORATIONS

ABUTMENTS

BioCam



BioCAM

STRAIGHT

[Hex]



23207031 1 mm
23207032 2 mm
23207033 3 mm

TCP



23607032 2 mm
23607033 3 mm

ANGLED



[Hex]

23231015 15° 2 mm
23231020 20° 2 mm

TITANIUM INTERFACE

[Hex]



23237905 0.5 mm
23237910 1 mm
23237920 2 mm
23237930 3 mm

TCP



23637920 2 mm
23637930 3 mm

SCREW-RETAINED RESTORATIONS

TAPERED



[Rotational]

23201031 10° 1 mm
23201032 10° 2 mm
23201033 10° 3 mm

ELEMENTS



13203313 Impression coping



41146000 Scanbody (TC2)



13203310 Healing cap



13205551 Analog



41005541 Digital analog



13209233 Titanium cylinder



13209133 Castable cylinder

TEMPORARY RESTORATIONS

ABUTMENTS

STRAIGHT



[Hex]

23232110 2 mm

TEMPORARY CYLINDERS



[Rotational]

13236110 Titanium
13235110 Castable

ANALOG



23205541

DIGITAL ANALOG



41235541

PROSTHETIC SCREWS



23232526 Hex 1.25 screw
232310151 For angled abutment
13202528 For tapered abutment

TCP



23632526 Hex H2-3 screw

MINI

STANDARD Platform restorations

HEALING ABUTMENTS

	23204810	1 mm
	23204820	2 mm
	23204830	3 mm
	23204840	4 mm
	23204850	5 mm
Ø 4	23204870	7 mm



Ø 5

23206036	3 mm
23206046	4 mm
23206056	5 mm
23206076	7 mm



Ø 6

AESTHETIC

23206030	3 mm
23206040	4 mm
23206050	5 mm
23206070	7 mm

SCANBODY



[Implant-level]
41236002 TC11

IMPRESSION COPINGS



[Open Tray]

23205503	Standard
23205506	Aesthetic
23255403	Long



[Open Tray]

23205504	Short
23205502	Long
23205505	Wax



[Close Tray]

23203724	Short
23203734	Long



[Close Tray]

23205600	Short
23205606	Long
23205603	Aesthetic

TCP

HEALING ABUTMENTS



23604820	2 mm
23604830	3 mm
23604840	4 mm
23604850	5 mm



IMPRESSION COPINGS

[Hex] [Open tray]
23605503



[Hex] [Close tray]
23603724 Corto
23603734 Largo



SCANBODY
[Implant-level]
41636002 TC11

SCREW-RETAINED RESTORATIONS

SINGLE



♥ticare
BioCAM

TITANIUM INTERFACE

[Hex]



23207905	0.5 mm
23207910	1 mm
23207920	2 mm
23207930	3 mm



TCP INTERFACE

23607920	2 mm
23607930	3 mm
23607940	4 mm

TITANIUM BASE FOR ZIRCONIUM

[Hex]



2320799* Titanium Base + CEREC
Scanbody (available in height
from 0 to 4 mm)
13207999 Scanbodies (3 units)

MULTIPLE

MULTIUNIT

STRAIGHT

[No hex]



23209902	2 mm
23209903	3 mm
23209904	4 mm
23209905	5 mm



HEALING SCREW
13201017

ANGLED

[Hex]



23201720	17°	2 mm
23201730	17°	3 mm
23201740	17°	4 mm
23203030	30°	3 mm
23203040	30°	4 mm
23203050	30°	5 mm



CYLINDERS

[Rot]

13209117	Titanium
13201751	Castable
13206710	Temporary titanium
13205710	Temporary castable

TITANIUM INTERFACE FOR ZIRCONIUM

13207917



SCANBODY
41146006 TC6



IMPRESSION COPING
13205217



ANALOG
13205117



DIGITAL ANALOG
41145117



♥ticare
BioCAM

CEMENT-RETAINED

ABUTMENT

BioCAM



STRAIGHT

[Hex]

23207000	0 mm Ø 4
23207010	1 mm
23207012	2 mm
23207013	3 mm
23207014	4 mm

[Hex]

23207613	3 mm Ø 5
23207614	4 mm

AESTHETIC

[Hex]

23207513	3 mm Ø 6
23207514	4 mm

23607012	2 mm	TCP
23607013	3 mm	
23607014	4 mm	

ANGLED

[Hex]

23201015	15° 1 mm Ø 4
23201315	15° 3 mm
23201515	15° 5 mm
23201020	20° 1 mm
23201320	20° 3 mm
23201520	20° 5 mm

[Hex]

23206015	15° 1 mm Ø 5
23206315	15° 3 mm
23206515	15° 5 mm
23206020	20° 1 mm
23206320	20° 3 mm
23206520	20° 5 mm

TITANIUM INTERFACE
FOR ZIRCONIUM

[Hex]

23207905	0.5 mm
23207910	1 mm
23207920	2 mm
23207930	3 mm

23607920	2 mm	TCP [Hex]
23607930	3 mm	
23607940	4 mm	

TITANIUM BASE FOR ZIRCONIUM

[Hex]

2320799*	Includes Scanbody CEREC (available in height from 0 to 4 mm)
13207999	Scanbodies (3 units)

REMOVABLE

ABUTMENT

BALL ABUTMENT

23208610	1 mm
23208620	2 mm
23208630	3 mm
23208640	4 mm
23208650	5.5 mm

* Plastic cap included. Ref 13208600

LOCATOR

23208681	1 mm
23208682	2 mm
23208683	3 mm
23208684	4 mm
23208685	5 mm

ELEMENTS

19008519	Processing package (2 units)
----------	---------------------------------

19008505	Impression coping Locator® (4units)
----------	--

DIGITAL ANALOG

41235501	
----------	--

ANALOG

23205501	Implant
23285501	Wrench
13205521	Ball abutment
19008530	Locator® (4 units)

TEMPORARY PROSTHESES

ABUTMENT

UCLA

[Hex] [Red]

23208151 Castable

[Non Hex] [White]

23207151 Castable

[Hex]

23204110 Titanium

[Non Hex]

23205110 Titanium

STRAIGHT

[Hex]

23202110 2 mm

[Non Hex]

23203110 2 mm

* Plastic cap included

PROSTHETIC SCREWS

DIRECT-TO-THE-IMPLANT

23202526	Hex 1.25 screw
23202508	Flat-connection screw
232010151	Angled abutment
23201726	Implant-angled Multiunit abutment retention screw

23602516	Hex. screw for TCP H2	TCP
23602526	Hex. screw for TCP H3-4	

NON DIRECT-TO-THE-IMPLANT

13202528	Tapered abutment-restoration prosthetic screw
13201727	Multiunit abutment-restoration prosthetic screw

LABORATORY

23282526	Lab screw (5 units)
----------	---------------------

STD

MAXI Platform restorations

HEALING ABUTMENTS



Ø 4

23205020	2 mm
23205030	3 mm
23205040	4 mm
23205050	5 mm
23205070	7 mm



Ø 5

23206536	3 mm
23206546	4 mm
23206556	5 mm
23206576	7 mm



Ø 6

AESTHETIC

23206530	3 mm
23206540	4 mm
23206550	5 mm
23206570	7 mm

SCANBODY



[Implant-level]
41236003 TC12

IMPRESSION COPINGS



[Open tray]

23205003	Standard
23205006	Long
23255403	Aesthetic



23205004	Short L=74mm
23205002	Long L=20mm
23205005	Wax L=29 mm



[Hex] [Close tray]

23205024	Short
23205034	Long



[No Hex] [Close tray]

23205056	Standard
23205066	Long
23205306	Aesthetic

TCP



HEALING SCREW

23605020	2 mm
23605030	3 mm
23605040	4 mm
23605050	5 mm



IMPRESSION COPING


[Hex] [Open tray]

23605003



[Hex] [Close tray]

23605024	Short
23605034	Long



SCANBODY

[Implant-level]

41636003 TC12

SINGLE



♥ticare
BioCAM

TITANIUM INTERFACE

[Hex]



23257905	0.5 mm
23257910	1 mm
23257920	2 mm
23257930	3 mm

TCP TITANIUM INTERFACE FOR ZIRCONIUM

[Hex]



23657920	2 mm
23657930	3 mm

TITANIUM BASE FOR ZIRCONIUM

[Hex]



2325799	Titanium Base + CEREC Scanbody (available in height from 0 to 4 mm)
13207999	Scanbodies (3 units)

MULTIPLE

MULTIUNIT

STRAIGHT

[No hex]



23259902	2 mm
23259903	3 mm
23259904	4 mm
23259905	5 mm



HEALING SCREW
13201017

ANGLED

[Hex]



23251720	17°	2 mm
23251730	17°	3 mm
23251740	17°	4 mm



23253030	30°	3 mm
23253040	30°	4 mm
23253050	30°	5 mm

CYLINDERS

[Rot]



13209117	Titanium
13201751	Castable
13206710	Temporary titanium
13205710	Temporary castable

TITANIUM INTERFACE FOR ZIRCONIUM

13207917



SCANBODY

41146006 TC06



IMPRESSION COPING

13205217



ANALOG

13205117



DIGITAL ANALOG

41145117



♥ticare
BioCAM

CEMENT-RETAINED RESTORATIONS

ABUTMENTS

BioCAM



STRAIGHT

[Hex]

23207050	0 mm Ø 4
23207051	1 mm
23207052	2 mm
23207053	3 mm
23207054	4 mm

[Hex]

23257613	3 mm Ø 5
23257614	4 mm

AESTHETIC

[Hex]

23257513	3 mm Ø 6
23257514	4 mm

TCP [Hex]

23607052	2 mm
23607053	3 mm
23607054	4 mm

ANGLED

[Hex]

23251015	15° 1 mm Ø 4
23251315	15° 3 mm
23251515	15° 5 mm
23251020	20° 1 mm
23251320	20° 3 mm
23251520	20° 5 mm

[Hex]

23256015	15° 1 mm Ø 5
23256315	15° 3 mm
23256515	15° 5 mm
23256020	20° 1 mm
23256320	20° 3 mm
23256520	20° 5 mm

TITANIUM INTERFACE

[Hex]

23257905	0.5 mm
23257910	1 mm
23257920	2 mm
23257930	3 mm

TCP [Hex]

23657920	2 mm
23657930	3 mm

TITANIUM BASE FOR ZIRCONIUM

[Hex]

2325799*	Includes Scanbody CEREC (available in height from 0 to 4 mm)
13207999	Scanbodies (3 units)

REMOVABLE

ABUTMENTS

BALL ABUTMENT



23258610	1 mm
23258620	2 mm
23258630	3 mm
23258640	4 mm
23258655	5.5 mm

* Retentive plastic cap. Ref 13208600

LOCATOR



23258681	0.73 mm
23258682	2 mm
23258683	3 mm
23258684	4 mm
23258685	5 mm

LOCATOR ELEMENTS



19008519 Processing pack (2 units)



19008505 Impression coping Locator® (4 units)

DIGITAL ANALOG



41235605

ANALOG



23205605 Implant



23285605 Wrench



13205525 Ball abutment



19008530 Locator® (4 units)

TEMPORARY PROTHESES & IMMEDIATE LOADING

ABUTMENTS

UCLA



[Hex] [Red]

23208155 Castable

[Non Hex] [White]

23207155 Castable



[Hex]

23254110 Titanium

[Non Hex]

23255110 Titanium

STRAIGHT



[Hex]

23252110 2 mm

[Non Hex]

23253110 2 mm

* Plastic cap included

PROSTHETIC SCREWS

DIRECT-TO-THE-IMPLANT



23205025 Hex 1.25 screw



23202509 Flat-connection screw



232510151 Angled abutment



23251726 Implant-angled multiunit abutment retention screw

TCP



23205025 Hex 1.25 screw

NON-DIRECT-TO-THE-IMPLANT



13202528 Tapered abutment-restoration prosthetic screw



13201727 Multiunit abutment-restoration prosthetic screw

LABORATORY



23285025 Lab screw (5 units)

MAXI



Spain

Teléfono: +34 983 211 312
sales@ticareimplants.com

MOZOGRAU S.A.

C/ Santiago López González, 7 · 47197 · Valladolid, Spain

www.ticareimplants.com
