




INHEX PROSTHETIC SYSTEM

♥ticare
i n H E X



GAP ZERO
TECHNOLOGY
ticare




Torque guide

		MINI	STD	MAXI
Cover screws		10 Ncm	10 Ncm	10 Ncm
Healing screws		10 Ncm	10 Ncm	10 Ncm
Impression copings		10 Ncm	10 Ncm	10 Ncm
Scanbodies		3 Ncm	3 Ncm	3 Ncm
Customized abutments		20 Ncm	30 Ncm	30 Ncm
Final abutments (direct to the implant)		20 Ncm	30 Ncm	30 Ncm
Components for 10° or 30° tapered abutments		15 Ncm	20 Ncm	20 Ncm
Components for 20° abutments (straight and angled multiunits)			15 Ncm	15 Ncm
Temporary abutments		10 Ncm	10 Ncm	10 Ncm
Screws for angled chimneys		Use the torque recommended by each manufacturer for the equivalent case of straight chimneys.		



Torque Ncm

MINI Platform PROSTHETIC COMPONENTS

HEALING ABUTMENTS




IMPRESSION COPINGS



TCP


HEALING ABUTMENT



IMPRESSION COPINGS



SCANBODY




SCREW-RETAINED RESTORATIONS


SINGLE




TITANIUM INTERFACE FOR ZIRCONIA



TCP TITANIUM INTERFACE FOR ZIRCONIA



SCANBODY




CEMENT-RETAINED RESTORATIONS

ABUTMENTS


BioCam




STRAIGHT





TCP



ANGLED




TITANIUM INTERFACE FOR ZIRCONIUM




SCREW-RETAINED RESTORATIONS

TAPERED ABUTMENTS



TCP




COMPONENTS




TEMPORARY RESTORATIONS

ABUTMENTS


STRAIGHT




ANALOG



DIGITAL ANALOG



PROSTHETIC SCREWS



MINI

STD Platform PROSTHETIC COMPONENTS

HEALING ABUTMENTS

	23204810 1 mm
	23204820 2 mm
	23204830 3 mm
	23204840 4 mm
	23204850 5 mm
Ø 4	23204870 7 mm



23206036 3 mm
23206046 4 mm
23206056 5 mm
23206076 7 mm



AESTHETIC

23206030 3 mm
23206040 4 mm
23206050 5 mm
23206070 7 mm

SCANBODY



[Implant-level]
41236002 TC11

IMPRESSION COPINGS



[Open Tray] [Non-hex]	23205503 Standard
	23205403 Long
	23205506 Aesthetic



[Open Tray]	23205504 Short
	23205502 Long
	23205505 Wax



[Close Tray] [Hex]	23203724 Short
	23203734 Long



[Close Tray] [Non-hex]	23205600 Short
	23205606 Long
	23205603 Aesthetic

TCP

HEALING ABUTMENTS



23604820 2 mm
23604830 3 mm
23604840 4 mm
23604850 5 mm



IMPRESSION COPINGS

[Hex] [Open tray]	23605503
-----------------------	----------



[Hex] [Close tray]	23603724 Corto
	23603734 Largo



SCANBODY
[Implant-level]
41636002 TC11

SCREW-RETAINED RESTORATIONS

SINGLE



TITANIUM INTERFACE FOR ZIRCONIA [Hex]

23207905 0.5 mm
23207910 1 mm
23207920 2 mm
23207930 3 mm



TCP	23607920 2 mm
	23607930 3 mm
	23607940 4 mm

TITANIUM INTERFACE FOR CEREC®

23207990 0 mm
23207991 1 mm
23207992 2 mm
23207993 3 mm
23207994 4 mm

MULTIPLE

TITANIUM INTERFACE FOR ZIRCONIA [Rot]

23207805 0.5 mm
23207810 1 mm
23207820 2 mm
23207830 3 mm



TCP	23607820 2 mm
	23607830 3 mm
	23607840 4 mm

MULTIPLE

MULTIUNIT

STRAIGHT

[No hex]



23209902 2 mm
23209903 3 mm
23209904 4 mm
23209905 5 mm



HEALING SCREW
13201017



TCP	23609902 2 mm
	23609903 3 mm
	23609904 4 mm
	23609905 5 mm

ANGLED

[Hex]



23201720 17° 2 mm
23201730 17° 3 mm
23201740 17° 4 mm



23203030 30° 3 mm
23203040 30° 4 mm
23203050 30° 5 mm



CYLINDERS

[Rot]

13209117 Titanium
13201751 Castable
13206710 Temporary titanium
13205710 Temporary castable

TITANIUM INTERFACE
FOR ZIRCONIUM
13207917



SCANBODY
41146006 TC6



IMPRESSION COPING
13205217



ANALOG
13205117



DIGITAL ANALOG
41145117



♥ticare
BioCAM

CEMENT-RETAINED

ABUTMENT


BioCAM



STRAIGHT




[Hex] Ø 4
23207000 0 mm
23207010 1 mm
23207012 2 mm
23207013 3 mm
23207014 4 mm




[Hex] Ø 5
23207613 3 mm
23207614 4 mm

AESTHETIC




[Hex] Ø 6
23207513 3 mm
23207514 4 mm



[Hex]
23607012 2 mm
23607013 3 mm
23607014 4 mm

ANGLED




[Hex] Ø 4
23201015 15° 1 mm
23201315 15° 3 mm
23201515 15° 5 mm
23201020 20° 1 mm
23201320 20° 3 mm
23201520 20° 5 mm




[Hex] Ø 5
23206015 15° 1 mm
23206315 15° 3 mm
23206515 15° 5 mm
23206020 20° 1 mm
23206320 20° 3 mm
23206520 20° 5 mm

TITANIUM INTERFACE FOR ZIRCONIA




[Hex]
23207905 0.5 mm
23207910 1 mm
23207920 2 mm
23207930 3 mm




[Non-Hex]
23207805 0.5 mm
23207810 1 mm
23207820 2 mm
23207830 3 mm

TCP



[Hex]
23607920 2 mm
23607930 3 mm
23607940 4 mm



[Non-Hex]
23607820 2 mm
23607830 3 mm
23607840 4 mm

TITANIUM INTERFACE
FOR CEREC®

23207990 0 mm
23207991 1 mm
23207992 2 mm
23207993 3 mm
23207994 4 mm

REMOVABLE

ABUTMENT

BALL ABUTMENT



23208610 1 mm
23208620 2 mm
23208630 3 mm
23208640 4 mm
23208650 5.5 mm



* Plastic cap included. Ref 13208600

LOCATOR



23208681 1 mm
23208682 2 mm
23208683 3 mm
23208684 4 mm
23208685 5 mm

ELEMENTS



19008519 Processing package (2 units)



190080282 Processing pack for bars (2 units)
--



19008505 Impression coping Locator® (4 units)



190085892 Locator for bar 2 mm (2 units)
--

DIGITAL ANALOG



41235501 Digital analog

ANALOGS



23205501 Implant analog



23285501 Wrench



13205521 Ball abutment



19008530 Locator® analog (4 units)

TEMPORARY ABUTMENTS

ABUTMENT

UCLA



[Hex] [Red]
23208151 Castable
[Non Hex] [White]
23207151 Castable



[Hex]
23204110 Titanium
[Non Hex]
23205110 Titanium

STRAIGHT



[Hex]
23202110 2 mm
[Non Hex]
23203110 2 mm

* Plastic cap included

PROSTHETIC SCREWS

DIRECT-TO-THE-IMPLANT



23202526 Hex 1.25 screw



23202508 Flat-connection screw



232010151 Angled abutment



23201726 Implant-angled Multiunit abutment retention screw
--



23602516 Hex. screw for TCP H2



23602526 Hex. screw for TCP H3-4

NON DIRECT-TO-THE-IMPLANT



13202528 Tapered abutment-restoration prosthetic screw
--



13201727 Multiunit abutment-restoration prosthetic screw
--

LABORATORY



23282526 Lab screw (5 units)

STD

MAXI Platform PROSTHETIC COMPONENTS

HEALING ABUTMENTS



23205020 2 mm
23205030 3 mm
23205040 4 mm
23205050 5 mm
23205070 7 mm



23206536 3 mm
23206546 4 mm
23206556 5 mm
23206576 7 mm



AESTHETIC

23206530 3 mm
23206540 4 mm
23206550 5 mm
23206570 7 mm

SCANBODY



[Implant-level]
41236003 TC12

IMPRESSION COPINGS



[Open tray]

23205003 Standard
23205007 Long
23255403 Aesthetic



23205004 Short L=14mm
23205002 Long L=20mm
23205005 Wax L=29 mm



[Hex] [Close tray]

23205024 Short
23205034 Long



[No Hex] [Close tray]

23205056 Standard
23205066 Long
23205306 Aesthetic

TCP



HEALING SCREW

23605020 2 mm
23605030 3 mm
23605040 4 mm
23605050 5 mm



IMPRESSION COPING


[Hex] [Open tray]

23605003



[Hex] [Close tray]

23605024 Short
23605034 Long



SCANBODY

[Implant-level]

41636003 TC12

SINGLE



TITANIUM INTERFACE FOR ZIRCONIA [Hex]

23257905 0.5 mm
23257910 1 mm
23257920 2 mm
23257930 3 mm



TCP

23657920 2 mm
23657930 3 mm
23657940 4 mm



TITANIUM INTERFACE FOR CEREC®

23257990 0 mm
23257991 1 mm
23257992 2 mm
23257993 3 mm
23257994 4 mm

MULTIPLE

TITANIUM INTERFACE FOR ZIRCONIA [Rot]

23257905 0.5 mm
23257910 1 mm
23257920 2 mm
23257930 3 mm



TCP

23657820 2 mm
23657830 3 mm
23657840 4 mm

MULTIPLE

MULTIUNIT

STRAIGHT

[No hex]



23259902 2 mm
23259903 3 mm
23259904 4 mm
23259905 5 mm



HEALING SCREW
13201017

ANGLED

[Hex]



23251720 17° 2 mm
23251730 17° 3 mm
23251740 17° 4 mm



23253030 30° 3 mm
23253040 30° 4 mm
23253050 30° 5 mm

CYLINDERS

[Rot]



13209117 Titanium
13201751 Castable
13206710 Temporary titanium
13205710 Temporary castable



TITANIUM INTERFACE FOR ZIRCONIUM
13207917



SCANBODY
41146006 TC6



DIGITAL ANALOG
13205217



ANALOG
13205117



IMPRESSION COPING
41145117

TCP



23659902 2 mm
23659903 3 mm
23659904 4 mm
23659905 5 mm



♥ticare
BioCAM

CEMENT-RETAINED RESTORATIONS

ABUTMENTS

BioCAM



STRAIGHT

[Hex] Ø 4

23207050 0 mm
23207051 1 mm
23207052 2 mm
23207053 3 mm
23207054 4 mm

[Hex] Ø 5

23257613 3 mm
23257614 4 mm

AESTHETIC

[Hex] Ø 6

23257513 3 mm
23257514 4 mm

TCP [Hex]

23607052 2 mm
23607053 3 mm
23607054 4 mm

ANGLED

[Hex] Ø 4

23251015 15° 1 mm
23251315 15° 3 mm
23251515 15° 5 mm
23251020 20° 1 mm
23251320 20° 3 mm
23251520 20° 5 mm

[Hex] Ø 5

23256015 15° 1 mm
23256315 15° 3 mm
23256515 15° 5 mm
23256020 20° 1 mm
23256320 20° 3 mm
23256520 20° 5 mm

TITANIUM INTERFACE FOR ZIRCONIA

[Hex]

23257905 0.5 mm
23257910 1 mm
23257920 2 mm
23257930 3 mm

[Non-Hex]

23257905 0.5 mm
23257910 1 mm
23257920 2 mm
23257930 3 mm

TCP

[Hex]

23657820 2 mm
23657830 3 mm
23657840 4 mm

[Non-Hex]

23657920 2 mm
23657930 3 mm

TITANIUM INTERFACE FOR CEREC®

23257990 0 mm
23257991 1 mm
23257992 2 mm
23257993 3 mm
23257994 4 mm

REMOVABLE

ABUTMENTS

BALL ABUTMENT

23258610 1 mm
23258620 2 mm
23258630 3 mm
23258640 4 mm
23258655 5.5 mm

* Retentive plastic cap. Ref 13208600

LOCATOR

23258681 0.73 mm
23258682 2 mm
23258683 3 mm
23258684 4 mm
23258685 5 mm

LOCATOR ELEMENTS

19008519 Processing package (2 units)
190080282 Processing pack for bars (2 units)
19008505 Impression coping Locator® (4 units)
190085892 Locator for bar 2 mm (2 units)

DIGITAL ANALOG

41235605 Digital analog

ANALOG

23205605 Implant analog
23285605 Wrench
13205521 Ball abutment
19008530 Locator® analog (4 units)

TEMPORARY ABUTMENTS & IMMEDIATE LOADING

ABUTMENTS

UCLA

[Hex] [Red]

23208155 Castable

[Non Hex] [White]

23207155 Castable

[Hex]

23254110 Titanium

[Non Hex]

23255110 Titanium

STRAIGHT

[Hex]

23252110 2 mm

[Non Hex]

23253110 2 mm

* Plastic cap included

PROSTHETIC SCREWS

DIRECT-TO-THE-IMPLANT

23205025 Hex 1.25 screw
23202509 Flat-connection screw
232510151 Angled abutment
23251726 Implant-angled multiunit abutment retention screw

23205025 Hex 1.25 screw TCP

NON-DIRECT-TO-THE-IMPLANT

13202528 Tapered abutment-restoration prosthetic screw
13201727 Multiunit abutment-restoration prosthetic screw

LABORATORY

23285025 Lab screw (5 units)

MAXI



Spain

Telephone: +34 983 211 312
sales@ticareimplants.com

MOZOGRAU S.A.

C/ Santiago López González, 7 · 47197 · Valladolid, Spain

www.ticareimplants.com
